Disposable Passive Accessory for Medtronic StealthStation® System

In the European Union, this device is not approved for direct contact with the central nervous system.

Intended Use

The Disposable Passive Accessory is intended to be used for anatomy palpation and registration in image guided surgeries with Medtronic StealthStation system. The device is sterile and designed for single use.

Contraindications

There are no known contraindications for this device.

Side Effects

Potential side effects associated with tumor resection procedures include: minor bleeding, hematoma without neurological compromise, or infection. Potential side effects associated with biopsy procedures include: minor bleeding, hematoma without neurological compromise, infection, or intraoperative brain bulge.

Adverse Events

Potential adverse events associated with tumor resection procedures, though not necessarily associated with use of this device, include but are not limited to: hemorrhage, hematoma with neurological compromise, seizure, brain swelling, infarction, disturbance of consciousness, motor deterioration, hemiparesis, aphasia, or other neurologic deficit.

Potential adverse events associated with biopsy procedures, though not necessarily associated with use of this device, include but are not limited to: fresh neurodeficit, pseudomeningocele, brain hemorrhage, or new onset seizure.

Warnings

- Prior to use, examine the device for damage, deterioration and inspect sterile barrier for breach. Do not use any compromised device. Abandon use of any device damaged during the procedure.
- This device is designed for single use and is supplied in sterile condition. Do not re-use or re-sterilize the product.
- During navigation, frequently confirm accuracy by touching the tip of the instrument on known anatomical points, including accuracy checkpoints, and comparing the position of the instrument tip in the image with its physical location.
- The device is MR UNSAFE.

Precautions

- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
- Do not re-use the devices
- Do re-sterilize the devices

Storage

Store in dry and clean conditions at ambient room temperature

Use

- Remove the pouch from the carton; carefully peel open the pouch and drop the sterile tray onto the sterile field or alternatively peel open the pouch and present the sterile tray to sterile personal.
- Carefully open the sterile tray, not to dislodge the loose Spherz markers.
- Verify the probe in the StealthStation interface.
 - Select the blunt probe from the instruments list in the procedure software.
 - Check that the probe has green status (it is being tracked by the StealthStation® system) on the screen.
 - Place the tip of the probe perpendicular in the reference frame divot.
 - Face the tracker array toward the camera.
 - Press the footswitch or wait for the system to verify.
- Use the probe to identify (localize) points on the patient anatomy or to simulate an instrument trajectory.

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IZI Medical Products, LLC
ILUMARK GmbH