IZI Medical

CE 0123

- EN**
3 **Duro-Ject™ Bone Cement Injector**
Instructions for Use
- DA**
4 **Duro-Ject™-knoglecementinjektor**
Brugsanvisning
- DE**
5 **Duro-Ject™ Knochenzementinjektor**
Gebrauchsanweisung
- EL**
6 **Εγχυτής συνδετικής ύλης οστών Duro-Ject™**
Οδηγίες Χρήσης
- ES**
8 **Inyector de cemento óseo Duro-Ject™**
Instrucciones de uso
- FR**
9 **Injecteur de ciment osseux Duro-Ject™**
Mode d'utilisation
- IT**
10 **Iniettore per cemento osseo Duro-Ject™**
Istruzioni per l'Uso
- NL**
11 **Duro-Ject™ Botcement-injector**
Gebruiksaanwijzing
- PT**
13 **Injector de cimento ósseo Duro-Ject™**
Instruções de Utilização
- SV**
14 **Duro-Ject™ Bencementsinjektor**
Bruksanvisning

T - D U R O - R E V 2

DURO-JECT™ BONE CEMENT INJECTOR

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Duro-Ject Bone Cement Injector is a manually powered 10 cc syringe with connecting tube used for placing bone cement into surgical sites.

INTENDED USE

The Duro-Ject Bone Cement Injector has been designed for placing bone cement into surgical sites.

CONTRAINDICATIONS

None known.

WARNINGS

- To avoid overfilling bone site and to stop the column of bone cement from advancing, user must turn plunger handle counterclockwise at least one full turn.
- If, during injection, you feel an increase in pressure and feel or see that no bone cement is being advanced out of the syringe, do not continue injection. Release the pressure by turning plunger handle counterclockwise at least one full turn and eliminate flow obstruction before continuing.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in bone access and infusion techniques.
- Follow all institutional protocols for the procedure.
- Review the Instructions for Use booklet prior to using this device.
- Injector is designed to inject viscous materials. Follow the bone cement manufacturer's instructions exactly.
- Turn the injector handle clockwise in a slow and controlled fashion.
- Monitor the entire cement injection under live fluoroscopic visualization.
- Due to the pressure that may be generated by this injector, safety glasses and/or a face shield are recommended.
- Do not over-pressurize the system; connecting tube may rupture.

INSTRUCTIONS FOR USE

1. Remove the injector/handle, syringe barrel, connecting tubes, funnels and extra o-ring assembly from the package.
2. Attach the connecting tube tightly to the Luer lock fitting on the tip of the syringe barrel.
3. Mix the bone cement according to the manufacturer's instructions.
4. Holding the syringe barrel with the tip down, place the funnel into the opening on the back of the barrel and pour the cement into the syringe barrel, being careful not to fill past the 10 ml mark.
5. While maintaining the vertical position of the syringe barrel, line up the plunger on the injector handle with the opening on the back of the syringe barrel.
6. Insert the plunger into the syringe barrel and turn the barrel clockwise until it locks into position.
7. Begin injecting cement by turning the plunger handle in a clockwise direction.

NOTE: One full turn of the handle delivers .20 cc of cement.
8. Advance cement through the connecting tube in small increments, until a small amount of cement is expelled, which will ensure the system is free of air.
9. Securely attach the connecting tube to the access needle.
10. Slowly turn the plunger handle, advancing the cement under live fluoroscopic visualization until an adequate fill has been achieved.
11. Stop the cement flow by turning the plunger handle in a counterclockwise direction at least one full turn.

NOTE: If it is necessary to reposition the needle within the bone, relieve pressure by turning the plunger handle counterclockwise at least one full turn before repositioning the needle.
12. Disconnect the connecting tube from the needle before removing the needle from the patient.

NOTE: If an additional injection is required, a second barrel, o-ring assembly, connecting tube, and funnel have been provided. Refer to steps 13 through 15 to complete additional injection.
13. Remove used o-ring assembly by turning plunger handle counterclockwise until o-ring assembly disengages from threaded plunger rod.
14. Place new o-ring assembly on the end of the threaded plunger rod, and push the o-ring assembly until it snaps in place. It might be necessary to rotate the plunger handle clockwise a few times before the o-ring assembly will fully engage.

15. Repeat steps 2 through 9 above to complete second injection.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local IZI Sales Representative for information on available literature

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IZI Medical Products, LLC
ILUMARK GmbH