



## **Instruction for Use**

Bone Drill  
Size: 11GA

Sterile  
For Single Use Only

## **Indication for Use**

The Bone Drill is intended for percutaneous access to the bone with the Vertebral Augmentation System.

## **Device Description**

The Bone Drill is packaged as a single-use, sterile, non-implantable device. It is a percutaneous surgical instrument designed to be used in conjunction with an appropriately sized Working Cannula, for the purpose of penetrating bone tissue for vertebral body access during Vertebral Augmentation procedures. The device consists of one primary component:

1. Bone Drill (Compatible with IZI Medical's 11GA Access Needle)

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

## **IMPORTANT:**

These instructions provide guidance to the experienced physician using the Bone Drill for accessing the vertebral body during a Vertebral Augmentation Procedure. This is not a reference to the general surgery technique. **Failure to properly follow the instructions may result in serious patient injury. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.**

## **Contraindications**

For the use of vertebral body access as determined by a licensed physician only. This product should be used by a physician familiar with the potential side effects of this procedure. Patients with bleeding disorders, or those receiving anti-coagulant therapies, should be evaluated for this procedure under physician judgment.

## **Precautions**

- Read the Instruction for Use prior to device operation.
- For single-patient use only.
- Product may present a biohazard and sharps hazard. Dispose of in accordance with applicable laws and regulations



**IZI Medical**

**Bone Drill Procedure:**

1. This instrument must be used in conjunction with an appropriately sized Working Cannula (not supplied) to provide a working channel to the vertebra.
2. Using a sterile technique, prepare the skin with antiseptic and drape.
3. Infiltrate the marked area with a local anesthetic, especially the periosteum.
4. Make a skin incision with a scalpel blade over the marked area.
5. Insert the 11GA access needle into the vertebral body using the Instruction for Use (IFU) for the access needle.
6. Insert the bone drill coaxially into the 11GA access needle.
7. Hold the Bone Drill with the proximal end in the palm and the index finger against the shaft near the tip. This position stabilizes the needle and allows better control.
8. Using gentle, but firm pressure, advance the Bone Drill into the bone. Rotate the needle in an alternating clockwise and counterclockwise motion while applying forward pressure. If the bone is very hard, lightly striking the handle with a mallet may advance the needle more easily.
9. Once inside the vertebral body, remove the Drill, leaving the Working Cannula in place within the pedicle.
10. Continue with the Vertebral Augmentation procedure per the associated Instructions for Use.

**Storage**

Store in a cool, dry place. Proper care should be taken to ensure that the product will not be damaged or rendered non-sterile.

**Limitation of Liability**

IN NO EVENT SHALL IZI MEDICAL BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, CONSEQUENTIAL, OR EXEMPLARY DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE PRODUCT, BASED UPON BREACH OF CONTRACT (INCLUDING BREACH OF WARRANTY).



**STERILE EO**

**Rx only**






**IZI Medical**

**Service Representatives and Requests for Information**

For service, technical support, requests for information, or reorder information, contact:



IZI Medical Products  
5 Easter Court, Suite J  
Owings Mills, MD 21117 USA  
+1 410 594-9403 (T)  
+1 410 594 0540 (F)  
[www.izimed.com](http://www.izimed.com)



ILUMARK GmbH  
Hohenlindner Str. 11c  
85622 Feldkirchen  
GERMANY

**Product Returns**

To return product, contact IZI Medical Inc. at (410) 594-9403