EN

Instruction for Use

OSTEO-SITE BONE FILLER INJECTION KIT
– Vertefix HV, 10GA / 11GA
Sterile
For Single Use Only

Indication for Use

Intended for the delivery of bone cement to bone with the Osteo-site Vertebral Balloon.

Device Description

The Bone Filler Injection Kit is packaged as single-use, sterile, non-implantable device. It is a percutaneous surgical instrument designed to be used in conjunction with an appropriately sized Kyphoplasty Working Cannula, for the purpose of penetrating bone tissue for vertebral body access during Kyphoplasty Vertebral Augmentation procedures. The device consists of the following components:

- 1. Bone Filler (Compatible with IZI Medical’s 10GA and 11GA Access Cannula, x6)
2. Osteo-force Injector
3. 2cc Syringes (x6)
4. Minimix Mixer CMM1
5. Vertefix HV Cement

Caution

Federal (USA) law restricts this device to sale by or on the order a physician (or properly licensed practitioner).

IMPORTANT:

These instructions provide guidance to the experienced physician using the Bone Filler Injection Kit for the delivery of bone cement to the vertebral body during a Vertebral Augmentation Procedure. This is not a reference to general surgery technique. Failure to properly follow the instructions may result in serious patient injury. The individual practitioner is responsible for the proper

procedure and techniques to be used with this device.

Contraindications

For the use of vertebral body access as determined by a licensed physician only. This product should be used by a physician familiar with the potential side effects of this procedure. Patients with bleeding disorders, or those receiving anti-coagulant therapies, should be evaluated for this procedure under physician judgment. If the flow of cement stops at any time during dispensing to determine and correct the cause the cause of blockage. Remove cement delivery system and bone filler device from the patient when to prevent them from becoming fixated within the vertebral body.

Warnings

None known

Precautions

- For single patient use only.
• Product may present a biohazard and sharps hazard. Dispose of in accordance with applicable laws and regulations.
• DO NOT use if package or component is damaged. DO NOT reuse, reprocess, or repackage this device. DO NOT modify any system component or accessory unless otherwise specified.
• Use only IZI-approved system components and accessories, unless otherwise specified.
• Only use with Polymethyl-Methacrylate (PMMA) based bone cement. DO NOT use with bone substitutes or bone void fillers. DO NOT overfill bone void.

Bone Filler Procedure

1. Follow the Instructions for Use (IFU), for the chosen Bone Access Tools to create access channel in the bone.
2. Follow the IFU for the chosen procedure to create a void or cavity in the vertebral body.
3. Begin the PMMA cement preparation process and follow the cement timings

IZI Medical Products
5 Easter Court, Suite J
Owings Mills, MD 21117 USA
+1 410 594-9403 (T)
+1 410 594 0540 (F)
www.izimed.com

ILUMARK GmbH
Hohenlindner Str. 11c
85622 Feldkirchen
Germany

stated on the cement IFU. Follow the images in Appendix for cement mixer instructions.

4. Note: Transfer cement from mixer to either syringes or Osteoforce injector. Osteoforce injector is preferred. Remove the Osteoforce plunger while filling it from the mixer.
5. Place the plunger back in the Osteoforce injector and turn the white knob to engage the threads of the Osteoforce injector.
6. Connect and fill bone filler from the Osteoforce injector.
7. Fill until completely full. Fill one at a time. Each turn of Osteoforce injector deploys 0.4cc of cement. Expect 2.5-3 turns to fill a bone filler (1cc).
8. Purge a small amount of cement to ensure no blockage. This will ensure the cement flows properly and safely. This also tests the connections of the system's components.
9. Allow time for the cement to reach a working state (doughy). Crank out a small amount of cement again to test the cement's viscosity to confirm the cement has reached a doughy state.
10. Using direct visualization or imaging guidance, advance the Bone Filler into the working cannula to the bone void. Verify placement of the cannula tip at the intended location using direct visualization or X-Ray imaging prior to delivery of bone cement.
11. The distal tip of the cannula has reached the end of the distal end of the cannula when the distal most marking on the cannula is inside the working cannula. Advance further if necessary.
12. Dispense slowly and consistently until the desired cement fill is achieved.
13. If the flow of cement stops at any time during dispensing, stop dispensing to determine and correct the cause of

blockage. Correct blockage or replace the Bone Filler.

14. When finished with cement fill, remove the system from the working cannula of the bone access tool. Tamp down with the Bone Filler, as needed, to verify the cement has set. Spin and rotate both the bone access tool and Bone Filler as required to prevent adhesion to the cement.
15. Note: If reusing the bone filler, clean the bone filler and inner stylet/pusher rod to ensure easy filling.
16. Continue with the Kyphoplasty Vertebral Augmentation procedure per associated Instructions for Use.

STORAGE

Store in a cool, dry place. Proper care should be taken to ensure that the product will not be damaged or rendered non-sterile.

References

These instructions for use are based on experience from physicians and (or) their published literature.

Limitation of Liability

IN NO EVENT SHALL IZI MEDICAL BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, CONSEQUENTIAL, OR EXEMPLARY DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE PRODUCT, BASED UPON BREACH OF CONTRACT (INCLUDING BREACH OF WARRANTY).

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