# SeeGrid®


Position indicator


- (EN) Instruction for use
- (SV) Bruksanvisning
- (NL) Handleiding
- (FR) Mode d'emploi
- (ES) Instrucciones de uso
- (IT) Istruzioni per l'uso
- (DE) Gebrauchsanweisung

**AprioMed**

(EN) Figure 1 (SV) Figur 1 (NL) Figuur 1 (DE) Abbildung 1 (FR) Figure 1  
(ES) Figura 1 (IT) Figura 1 (DE) Abbildung 1

(EN) Figure 2 (SV) Figur 2 (NL) Figuur 2 (DE) Abbildung 2 (FR) Figure 2  
(ES) Figura 2 (IT) Figura 2 (DE) Abbildung 2

(EN) Figure 3 (SV) Figur 3 (NL) Figuur 3 (DE) Abbildung 3 (FR) Figure 3  
(ES) Figura 3 (IT) Figura 3 (DE) Abbildung 3

# ENGLISH

## DEVICE DESCRIPTION

The system consists of the following parts:

1 SeeGrid®, packaged in a sterile bag.

## INDICATION FOR USE

To provide well defined position indicators following the body contour.

## WARNINGS (W)

⚠ **W1.** When using SeeGrid® together with other devices, please check compatibility of the tool with SeeGrid®.

⚠ **W2.** Please note that SeeGrid® adds 2 mm thickness to the skin surface.

⚠ **W3.** Make sure that both ends of SeeGrid® are visible along the body contour.

## PRECAUTIONS (P)

⚠ **P1.** SeeGrid® should only be used by, or under the supervision of, a physician trained in the appropriate procedures.

⚠ **P2.** Carefully inspect SeeGrid® prior to use to verify that it is undamaged. Check that the sterile packaging is intact.

⚠ **P3.** The system is designed for single patient use only. Do not re-use or re-sterilize. Store in a dry, cool, dark place.

⚠ **P4.** After use, the grid may be a potential biohazard. Handle and dispose of SeeGrid® in accordance with accepted medical practise and applicable local and national laws and regulations.

⚠ **P5.** SeeGrid® is not validated for use on breached or compromised skin.

## DIRECTIONS FOR USE

1. Place SeeGrid® on the skin surface of the body. The flexible material makes it easy to follow any body contour (see figure 1).

2. Use tape to fasten SeeGrid® on the skin. Make sure that it is securely fixed before proceeding (see figure 2).

3. A control scan must be taken in order to verify that all indicators are visible on the picture (see figure 3).

## OTHER INFORMATION

Serious incident that has occurred in relation to the device shall be reported in accordance with national laws and regulations. Serious incidents shall be reported to AprioMed AB (legal manufacturer) and to the appropriate regulatory agency such as the competent authority of the Member State (EU/EEA) in which the user/patient is established.

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SeeGrid® is a registered trademark of AprioMed AB.
Mosaic Surgical Ltd