
Gangi-SoftGuard[®]

Coaxial Biopsy Needle


- (EN) Instruction for use
- (SV) Bruksanvisning
- (NL) Handleiding
- (FR) Mode d'emploi
- (ES) Instrucciones de uso
- (IT) Istruzioni per l'uso
- (DE) Gebrauchsanweisung
- (NO) Bruksanvisning


AprioMed

(EN) Figures (SV) Figurer (NL) Afbeeldingen (FR) Figure (ES) Figuras (IT) Figura (DE) Abbildungen (NO) Figurer

1

2

(EN) Unlock

(SE) Lås upp

(NL) Ontgrendelen

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(ES) Desbloquee

(IT) Sbloccare

(DE) Öffnen

(NO) Lås opp

(EN) Remove

(SE) Dra ut

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(FR) Retirer

(ES) Retire

(IT) Rimuovere

(DE) Entfernen

(NO) Fjern

(EN) Insert

(SE) Sätt in

(EL) Inbrengein

(FR) Insérer

(ES) Inserte

(IT) Inserire

(DE) Einführen

(NO) Sett inn

(EN) Lock

(SE) Lås

(NL) Vergrendel

(FR) Verrouiller


(ES) Bloqueo

(IT) Bloccare


(DE) Schließen

(NO) Lås


3

4

5

6

(EN) Gangi-SoftGuard® consists of the following parts: (SE) Gangi-SoftGuard® består av: (NL) De Gangi-SoftGuard® bestaat uit de volgende onderdelen: (FR) L'aiguille Gangi-SoftGuard® est constituée des éléments suivants: (ES) Gangi-SoftGuard® consta de las siguientes partes: (IT) Gangi-SoftGuard® è composto dalle seguenti parti: (DE) Die Gangi-SoftGuard® besteht aus den folgenden Teilen: (NO) Gangi-SoftGuard® består av følgende deler:

(EN) Blunt Stylet (SE) Trubbig Mandrång (NL) Stompe stylet (FR) Mandrin arrondi (ES) Estilete Punta Roma (IT) Stiletto smussato (DE) Stillett mit abgestumpfter Spitze (NO) Butt stillett


(EN) Trocar Stylet (SE) Trokarmandrång (NL) Trocar stylet (FR) Mandrin pointu (ES) Estilete Punta Trocar (IT) Stiletto Trocar (DE) Stillett mit Trokar-Spitze (NO) Trocar stillett

(EN) Needle (SE) Nål (NL) Naald (FR) Aiguille (ES) Aguja (IT) Ago (DE) Nadel (NO) Nål

(EN) Reference Guide (SE) Referensguide (NL) Referentie gids (FR) Guide de référence (ES) Guía de Referencia (IT) Guida di riferimento (DE) Bedienungsanleitung (NO) Referansemal

(EN) Table 1. (SV) Tabell 1. (NL) Tabel 1. (FR) Tableau 1. (ES) Tabla 1. (IT) Tabella 1. (DE) Tabelle 1. (NO) Tabell 1.

REF	Gauge x Length Gauge x Längd Gauge x Lengte Gauge x Long Calibre x Longitud Gauge x Lunghezza Gauge x Länge Gauge x Lengde	Needle ID ID Nål Naald Binnen Ø Diam Int Aiguille DI Aguja Diam Int Ago Nadel Innen Ø Nål ID	Suitable biopsy device size (Gauge/Length) Storlek på passande biopsiinstrument (Gauge x Längd) Geschikte maat biopsie naald (Gauge x Lengte) Taille du dispositif de biopsie approprié (Gauge x Long) Tamaño del dispositivo de Biopsia (Calibre x Longitud) Misura idonea del dispositivo di biopsia (Gauge x Lunghezza) Geeignete Größe des Biopsie-Produktes (Gauge x Länge) Egnet størrelse biopsienhet (Gauge x Lengde)	≥ 14.1 cm ≥ 9.1 cm ≥ 14.1 cm ≥ 19.1 cm ≥ 9.1 cm ≥ 14.1 cm ≥ 19.1 cm ≥ 9.1 cm ≥ 14.1 cm ≥ 19.1 cm
● SG12-121	12G x 12.1 cm	2.48 mm	13G	≥ 14.1 cm
● SG13-071	13G x 7.1 cm	2.18 mm	14G	≥ 9.1 cm
● SG13-121	13G x 12.1 cm	2.18 mm	14G	≥ 14.1 cm
● SG13-171	13G x 17.1 cm	2.18 mm	14G	≥ 19.1 cm
● SG15-071	15G x 7.1 cm	1.73 mm	16G	≥ 9.1 cm
● SG15-121	15G x 12.1 cm	1.73 mm	16G	≥ 14.1 cm
● SG15-171	15G x 17.1 cm	1.73 mm	16G	≥ 19.1 cm
● SG17-071	17G x 7.1 cm	1.33 mm	18G	≥ 9.1 cm
● SG17-121	17G x 12.1 cm	1.33 mm	18G	≥ 14.1 cm
● SG17-171	17G x 17.1 cm	1.33 mm	18G	≥ 19.1 cm

ENGLISH

INDICATION FOR USE

As a guiding needle in obtaining core biopsy samples from soft tissue such as liver, kidney, spleen, lymph nodes and various soft tissue lesions.

CONTRAINDICATIONS

Not intended for use in bone.

WARNINGS (W)

- ⚠ **W1.** Good medical judgement should be exercised in considering biopsy on patients who are receiving anticoagulant therapy or who have bleeding disorders. Use in these patient populations may result in severe bleeding.
- ⚠ **W2.** In case of bending of any part of the coaxial system, stop the procedure and retract the bent part simultaneously with the rest of the device. Removing bent parts only may cause breakage and result in materials remaining in the body.
- ⚠ **W3.** Gangi-SoftGuard® is intended for use in image-guided procedures only. Verify the device location using imaging techniques after every device adjustment. Do not use Gangi-SoftGuard® in Magnetic Resonance Imaging (MRI) equipment. The device is magnetic and exposure to MRI may cause unintended device movement.
- ⚠ **W4.** One potential complication of coaxial image-guided biopsy is inadvertent perforation of adjacent structures. Inadvertent perforation may result in site specific complications which may consist of bleeding, including bleeding from the large artery, hematoma, hemorrhage, infection, pain, cardiac tamponade, pneumothorax, hemoptysis and hemothorax.

PRECAUTIONS (P)

- ⚠ **P1.** Before use together with other equipment, check the compatibility in length and gauge together with Gangi SoftGuard® using the reference guide included in the package.
- ⚠ **P2.** Before use, inspect the device for damaged point, bent shaft or other imperfections that would prevent proper function. If the needle components are damaged or bent or do not move together easily, DO NOT USE.
- ⚠ **P3.** After use, the components may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.
- ⚠ **P4.** Gangi SoftGuard® must be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of core needle biopsy, in particular, those relating to the specific organ being biopsied.

DIRECTION FOR USE

Preparation

Before use together with other biopsy instruments check the compatibility in length and gauge together with Gangi SoftGuard® Coaxial Needle or by using the Reference Guide provided. Confirm that the distal end of the biopsy instrument extends through the Gangi SoftGuard® Coaxial Needle tip. Compatible biopsy device gauge sizes and needle length are described in Table 1. Use aseptic technique while removing Gangi SoftGuard from the package. Before use, inspect the device for damaged point, bent shaft or other imperfections that would prevent proper function. If the needle components are damaged or bent or do not move together easily, DO NOT USE. ⚠ **W1, W2**

Biopsy procedure

1. Determine which stylet is preferred (Figure 1). If the blunt stylet is preferred, exchange stylets (Figure 2). The blunt tip stylet is used when advancing through soft tissue and around vasculature or other organs and reduce the risk of damage. Function of the blunt stylet (Figure 3).
2. The biopsy procedure must be performed using appropriate aseptic technique. Insert device superficially (Figure 4). Monitor position using imaging technique. Note: If trocar stylet is not used, make a small skin incision using a scalpel to make penetration of skin easier. ⚠ **W3**
3. If trocar stylet is used, exchange to blunt stylet. Advance needle proximal to lesion (Figure. 5). Monitor position using centimeter markings on the needle and imaging technique. Note: If necessary, stylets can be exchanged during procedure (Figure 2). ⚠ **W2, W3, W4**
4. Remove stylet. Perform biopsy as per normal routine (Figure 6). ⚠ **W3**
5. For multiple samples, reinsert stylet and reposition Gangi-SoftGuard using image guidance. Repeat step 1, 2, 3 and 4. ⚠ **W2, W3** Recommendation: When collecting multiple samples, wipe the core tissue biopsy needle with sterile moist gauze prior to reinsertion into the cannula needle. This will aid the proper movement of the core tissue biopsy needle within the cannula needle.
6. When the procedure is completed, Reinsert the stylet in the needle and remove the device.

OTHER INFORMATION

Serious incident that has occurred in relation to the device shall be reported in accordance with national laws and regulations. Serious incidents shall be reported to AprioMed AB (legal manufacturer) and to the appropriate regulatory agency such as the competent authority of the Member State (EU/EEA) in which the user/patient is established.

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Rx Only

(EN) Caution: Federal (U.S.A) law restricts this device to sale by or on the order of a physician.

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(EN) Single use. Do not re-use. After use, the product performance and characteristics may degrade, and can cause infection in patient.

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(EN) Do not re-sterilize. Re-sterilization may have impact on product characteristics and performance.

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(EN) Do not use if package is damaged. Device is delivered sterile. Product in damaged package may be contaminated and can cause infection in patient.

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(EN) MR unsafe. Do not use AprioCore plus in Magnetic Resonance Imaging (MRI) equipment. The device is magnetic and exposure to MRI may cause unintended device movement.

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Gangi-SoftGuard® is a registered trademark of AprioMed AB.
Mosaic Surgical Ltd