- EN**
2
Kopans Breast Lesion Localization Needle
Instructions for Use
- DA**
2
Kopans lokaliseringskanyle til læsioner i brystet
Brugsanvisning
- DE**
3
Kopans Brustläsions-Lokalisationskanüle
Gebrauchsanweisung
- EL**
4
Βελόνα εντοπισμού βλαβών μαστού Kopans
Οδηγίες χρήσης
- ES**
5
Aguja de localización de lesiones mamarias Kopans
Instrucciones de uso
- FR**
6
Aiguille de Kopans pour localisation de lésion mammaire
Mode d'emploi
- IT**
7
Ago per la localizzazione delle lesioni mammarie Kopans
Istruzioni per l'uso
- NL**
8
Kopans lokalisatienaald voor borstlaesies
Gebruiksaanwijzing
- PL**
9
Igła typu Kopans do lokalizacji zmian w piersi
Instrukcja użycia
- PT**
10
Agulha para localização de lesões mamárias Kopans
Instruções de utilização
- SV**
11
Kopans nål för bröstlesionslokalisering
Bruksanvisning
- ZH**
12
Kopans乳房病灶定位针
使用说明

KOPANS BREAST LESION LOCALIZATION NEEDLES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Kopans Breast Lesion Localization Needle consists of a bevel-tip needle with centimeter markings and a hookwire. The hookwire features a burnish mark that allows the operator to ensure that the hookwire is within the needle tip during needle manipulation.

INTENDED USE

The Kopans Breast Lesion Localization Needle is intended for preoperative marking of nonpalpable breast lesions.

CONTRAINDICATIONS

None known

WARNINGS

- Under no circumstances should a hookwire engaged in tissue be pulled out without surgical intervention.

PRECAUTIONS

- This product is made of ferrous material and has not been tested for safety in the MR environment.
- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for localization of breast lesions should be employed.
- Following placement of the hookwire, the portion protruding outside of the breast should be bent and taped to the skin to prevent inadvertent movement.
- Final hookwire position should be confirmed by appropriate imaging modality.
- The hookwire should be used as a guide for the surgeon, not a retractor.

POTENTIAL ADVERSE EVENTS

- Infection
- Internal bleeding and/or bleeding at the puncture site
- Pneumothorax
- Puncture or injury to target organ or to a nearby organ that is traversed by the needle.
- Hemoptysis
- Peritonitis
- Damage of tissue adjacent to the needle

INSTRUCTIONS FOR USE

NOTE: The burnish mark on the hookwire should be seen just outside the needle hub to ensure that the hookwire is within the needle tip during needle manipulation.

1. Introduce the needle into the lesion.
2. Check the position of the needle.
3. Advance and release the hookwire.
4. Remove the needle.
5. Bend the hookwire protruding from the breast and tape flat to the skin.
6. Verify final hookwire position.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© IZI Medical Products 2019
Owings Mills, MD 21117 U.S.A.
ILUMARK GmbH