Bonopty[®]

Bone Biopsy System



- (EN) Instruction for use
- (SV) Bruksanvisning
- (NL) Handleiding
- (FR) Mode d'emploi
- (ES) Instrucciones de uso
- (IT) Istruzioni per l'uso
- (DE) Gebrauchsanweisung
- (NO) Bruksanvisning


AprioMed

(EN) Contents (SV) Innehåll (NL) Inhoud (FR) Contenu (ES) Contenido (IT) Riepilogo (DE) Inhalt (NO) Innhold

(EN) TABLE 1A. PENETRATION SET (SV) TABELL 1A. PENETRATIONSSET (NL) TABEL 1A. PENETRATIESET (FR) TABLEAU 1A. MATÉRIEL DE PONCTION (ES) TABLA 1A. CONJUNTO DE PENETRACIÓN (IT) TABELLA 1A. SET DI PENETRAZIONE (DE) TABELLE 1A. PENETRATIONSSSET (NO) TABELL 1A. PENETRASJONSSET

REF	10-1062	10-1072	10-1082
(EN) Penetration Cannula with stylet (SV) Penetrationsrör med mandrin (NL) Penetratiecanule met stilet (FR) Canule de ponction avec mandrin (ES) Cánula de penetración con estilete (IT) Cannula di penetrazione con mandrino (DE) Penetrationskanüle mit Mandrin (NO) Penetrasjonsnål med stilet			
(EN) OD (SV) OD (NL) OD (FR) DE (ES) DE (IT) Diam. est. (DE) OD (NO) OD	2.10 mm		
(EN) ID (SV) ID (NL) ID (FR) DI (ES) DI (IT) Diam. int. (DE) ID (NO) ID	1.80 mm		
(EN) Length (SV) Längd (NL) Lengte (FR) Longueur (ES) Longitud (IT) Lunghezza (DE) Länge (NO) Lengde	L = 65 mm	L = 95 mm	L = 130 mm
(EN) Drill (SV) Borr (NL) Boor (FR) Foret (ES) Taladro (IT) Punta (DE) Bohrer (NO) Bor			
(EN) OD (SV) OD (NL) OD (FR) DE (ES) DE (IT) Diam. est. (DE) OD (NO) OD	1.70 mm		
(EN) Length (SV) Längd (NL) Lengte (FR) Longueur (ES) Longitud (IT) Lunghezza (DE) Länge (NO) Lengde	L = 92 mm	L = 122 mm	L = 160 mm
(EN) Depth Gauge (SV) Djupanslag (NL) Dieptemeter (FR) Butée de profondeur (ES) Calibrador de profundidad (IT) Calibro di profondità (DE) Tiefenmessgerät (NO) Dybdemåler			

(EN) TABLE 1B. BIOPSY SET (SV) TABELL 1B. BIOPSISET (NL) TABEL 1B. BIOPSIESET (FR) TABLEAU 1B. MATÉRIEL DE BIOPSIE (ES) TABLA 1B. CONJUNTO DE BIOPSIA (IT) TABELLA 1B. SET PER BIOPSIA (DE) TABELLE 1B. BIOPSIE-SET (NO) TABELL 1B. BIOP-SISETT

REF	10-1063	10-1073	10-1083
(EN) Biopsy Cannula with stylet (SV) Biopsirör med mandrin (NL) Biopsiecanule met stilet (FR) Canule pour biopsie avec mandrin (ES) Cánula de biopsia con estilete (IT) Cannula per biopsia con mandrino (DE) Biopsiekanüle mit Mandrin (NO) Biopsikanyle med stilet			
(EN) OD (SV) OD (NL) OD (FR) DE (ES) DE (IT) Diam. est. (DE) OD (NO) OD	1.70 mm		
(EN) ID (SV) ID (NL) ID (FR) DI (ES) DI (IT) Diam. int. (DE) ID (NO) ID	1.30 mm		
(EN) Length (SV) Längd (NL) Lengte (FR) Longueur (ES) Longitud (IT) Lunghezza (DE) Länge (NO) Lengde	L = 130 mm	L = 160 mm	L = 195 mm
(EN) Depth Gauge (SV) Djupanslag (NL) Dieptemeter (FR) Butée de profondeur (ES) Calibrador de profundidad (IT) Calibro di profondità (DE) Tiefenmessgerät (NO) Dybdemåler			
(EN) Ejector Pin (SV) Utstötartift (NL) Uitwerpen (FR) Éjecteur de carotte (ES) Pasador del eyector (IT) Eiettore (DE) Ausstoßer-Pin (NO) Utløsningspinne			


(EN) TABLE 1C. EXTENDED DRILL (SV) TABELL 1C. FÖRLÄNGD BORR (NL) TABEL 1C. VERLENGDE BOOR (FR) TABLEAU 1C. FORET ALLONGÉ (ES) TABLA 1C. TALADRO AMPLIADO (IT) TABELLA 1C. PUNTA ESTESA (DE) TABELLE 1C. LANGBOHRER (NO) TABELL 1C. UTVIDET BOR

REF	10-1074
(EN) Drill (SV) Borr (NL) Boor (FR) Foret (ES) Taladro (IT) Punta (DE) Bohrer (NO) Bor	
(EN) OD (SV) OD (NL) OD (FR) DE (ES) DE (IT) Diam. est. (DE) OD (NO) OD	1.70 mm
(EN) Length (SV) Längd (NL) Lengte (FR) Longueur (ES) Longitud (IT) Lunghezza (DE) Länge (NO) Lengde	L = 160 mm

(EN) Table 2. Protuding lenght (EN) Tabell 2. Utskjutande längd (EN) Tabel 2. Uitstekende lengte (FR) Tableau 2. Longueur de dépassement de la partie aiguisée (ES) Tabla 2. Longitud que sobresale (IT) Tabella 2. Lunghezza sporgente (DE) Tabelle 2. Überstehende Länge (NO) Tabell 2. Utstikkende lengde

	(EN) Protruding length of the biopsy cannula (SV) Utskjutande längd på biopsiröret (EN) Uitstekende lengte van de biopsiecanule (FR) Longueur de dépassement de la partie aiguisée de la canule de biopsie (ES) Longitud de la cánula de biopsia que sobresale (IT) Lunghezza sporgente della cannula per biopsia (DE) Überstehende Länge der Biopsiekanüle (NO) Utstikkende lengde på biopsikanylen			(EN) Protruding length of the extended drill (SV) Utskjutande längd på den förlängda borren (EN) Uitstekende lengte van de verlengde boor (FR) Longueur de dépassement du foret allongé (ES) Longitud del taladro ampliado que sobresale (IT) Lunghezza sporgente della punta estesa (DE) Überstehende Länge des Langbohrers (NO) Utstikkende lengde på det utvidede boret
	10-1063 Biopsy Set 13 cm	10-1073 Biopsy Set 16 cm	10-1083 Biopsy Set 19.5 cm	10-1074 Extended drill 16 cm
10-1062 Penetration Set 6.5 cm	X = 45 mm	X = 75 mm	X = 115 mm	X = 75 mm
10-1072 Penetration Set 9.5 cm	X = 10 mm	X = 45 mm	X = 85 mm	X = 45 mm
10-1082 Penetration Set 13 cm	--	X = 10 mm	X = 45 mm	--

(EN) Figure 1. (SV) Figur 1. (NL) Figuur 1. (FR) Figure 1. (ES) Figura 1. (IT) Figura 1. (DE) Abbildung 1. (NO) Figur 1.

ENGLISH

DEVICE DESCRIPTION

See Table 1A-1C and 2.

INDICATION FOR USE

Biopsy of bone lesions.

CONTRAINDICATIONS

The Bonoptoy® Coaxial Bone Biopsy System with Eccentric Drill is contraindicated for use on patients who are receiving heavy anticoagulant therapy or who have a severe bleeding problem.

WARNINGS (W):

⚠ **W1.** When using devices from other manufacturers, first check the compatibility (length and gauge) of the tool with the Bonoptoy® components.

⚠ **W2.** Use finger force only during the biopsy procedure. Do not use extra tools to manipulate the components of the Bonoptoy® System.

⚠ **W3.** Not suitable for very sclerotic or purely lytic lesions. If a target or trajectory proves to be too hard to sample or penetrate using the Bonoptoy® Biopsy Cannula, stop the procedure and remove the Cannula. It may be possible to sample using the Bonoptoy® Extended Drill instead; see section Bonoptoy® Extended Drill.

⚠ **W4.** The Penetration Cannula is designed for anchorage, and is not intended for passage through thick intact cortical bone.

⚠ **W5.** When the procedure is completed, do not remove the cannula without the support of the stylet or drill. Otherwise there is a risk that the cannula may break and part of it remains in the bone.

⚠ **W6.** Bonoptoy is intended for use in image-guided procedures only. Verify the device location using imaging techniques after every device adjustment. Do not use Bonoptoy in Magnetic Resonance Imaging (MRI) equipment. The device is magnetic and exposure to MRI may cause unintended device movement.

PRECAUTIONS (P):

⚠ **P1.** The Bonoptoy® Coaxial Bone Biopsy System with Eccentric Drill should be used only by or under the supervision of a physician trained in biopsy or radiological interventions.

⚠ **P2.** Carefully inspect the contents of the system prior to use to verify that all parts are present and undamaged.

⚠ **P3.** The Bonoptoy® System is designed for single patient use only. After one use, components of the system may become dull or damaged, making a biopsy more difficult and less successful.

⚠ **P4.** Be careful not to bend the components of the Bonoptoy® System during manipulation.

⚠ **P5.** Contaminated products should be disposed utilizing standard hospital procedures for biohazardous waste.

⚠ **P6.** During drilling, do not push the Cannula forward over the Drill, this hinder the "wobbling" movement of the eccentric Drill and prevent the transportation of bone material out from the hole

DIRECTION FOR USE

1. The BONOPTY® Penetration Set ⚠ **W1 - W6**, ⚠ **P1- P6**

The following is a method for penetrating cortical bone and creating a fixed pathway for coaxial sampling using the BONOPTY® Penetration Set.

a. Apply local anaesthetic.

b. Make a small skin incision prior to the insertion of the Penetration Cannula and Stylet.

c. Insert the Penetration Cannula and Stylet through the soft tissue towards the bone surface. (Figure 1.1)

Note: If necessary, the stylet can be removed to allow additional anaesthetic injections through the opening of the Penetration Cannula. (Figure 1.2)

d. When the bone is reached, exchange the stylet for the Drill (Figure 1.3)

Note: During this exchange, the Penetration Cannula must be firmly held against the bone so that a stable position is maintained.

e. Rotate the Drill clockwise a few turns using moderate axial pressure until the tip of the Drill is anchored in the bone. (Figure 1.4)

f. Verify the position and direction of the Drill and Penetration Cannula using diagnostic imaging. ⚠ **W6**

Note: When necessary, it is now easy to correct the position and angulation of the Drill and Penetration Cannula.

Note: The tip of the Drill protrudes a maximum of 8mm beyond the tip of the Penetration Cannula.

Note: During drilling, do not push the Cannula forward over the Drill, as this may hinder the "wobbling" movement of the eccentric Drill and prevent the transportation of bone material out from the hole. (Figure 1.4)

g. Rotate the Drill clockwise further into the bone to the desired depth just beyond the cortical bone. (Figure 1.5)

h. Advance the Penetration Cannula over the Drill by careful clockwise rotation and anchor the Penetration Cannula 2-3 millimetres into the bone. Δ **W4**. The Penetration Cannula is designed for anchorage, and is not intended for passage through thick intact cortical bone

i. Verify the position and direction of the Penetration Cannula and the Drill using diagnostic imaging.

j. Withdraw the Drill a few centimetres and verify the anchorage of the Penetration Cannula by careful manipulation.

k. Remove the Drill completely. The anchored Penetration Cannula may now be used as a fixed pathway for biopsy sampling. (Figure 1.6)

2. The BONOPTY® Biopsy Set Δ **W1 - W6**, Δ **P1- P6**

The following is a method for sampling bone lesions using the BONOPTY Biopsy Set through the fixed pathway provided by the BONOPTY Penetration Set. Not suitable for very sclerotic or purely lytic lesions.

a. The Biopsy Cannula protrudes a maximum of x cm beyond the Penetration Cannula tip (x= see Figure B). The protruding length can be shortened up to 4 cm by using the Depth Gauge (break indicators every 5mm):

- Fit the Depth Gauge onto the proximal opening on the hub of the Penetration Cannula. (Figure 1.7)

- Break the Depth Gauge at a chosen letter. (Figure 1.8)

b. Insert the Biopsy Cannula with its stylet through the Penetration Cannula. Rotate the Biopsy Cannula clockwise with moderate force until the biopsy target is reached. Break off the number of indicators corresponding to the desired length of the sample. (Figure 1.9 - 1.10)

Note: The stylet of the Biopsy Set should not be withdrawn until the target has been reached to avoid filling the Biopsy Cannula with drilled bone fragments or non-target tissue on the way to the target.

If the pathway to the target is difficult to penetrate using the Biopsy Cannula with stylet, it is possible to switch to the Extended Drill in order to gain access to the target.

c. Withdraw the stylet and obtain a sample by rotating the Biopsy Cannula. (Figure 1.11) of the stylet or drill. Otherwise there is a risk that the cannula may break and part of it remains in the bone. Δ **W5**

Note: The Biopsy Cannula seldom obtains samples longer than 20mm. If more specimen material is required, it is advisable to perform repeated coaxial biopsies at different depths in large lesions. Warning: If a target proves to be too hard to sample using the Biopsy Cannula, stop the procedure and remove the Biopsy Cannula. It may be possible to sample using the

Extended Drill instead; see section 3. Δ **W3**

d. Rotate the Biopsy Cannula while withdrawing it. (Figure 1.12)

e. Because of the inward reduction of the Biopsy Cannula tip, the specimen must be ejected from the hub end in order to avoid crush artefacts. Insert the Ejector Pin into the Cannula tip forcing the sample out through the proximal opening at the Cannula hub. (Figure 1.13)

Note: Be careful not to damage the sharp cutting edge of the Cannula when inserting the Ejector Pin.

f. When the procedure is completed, do not remove the penetration cannula without the support

3. The BONOPTY® Extended Drill Δ **W1-W2**, Δ **W4-W6**, Δ **P1-P6**

The following is a method for sampling very sclerotic lesions using the BONOPTY® Extended Drill.

Note: The Drill may also be used through the Penetration Cannula to reach a target if the pathway is difficult to penetrate using the Biopsy Cannula with stylet; see section 2.

a. Insert the Extended Drill into the Penetration Cannula. Drill protrudes a maximum of x cm beyond the Penetration Cannula tip (x= see Table 2).

b. Turn the Extended Drill a few clockwise rotations to fill the helices of the Extended Drill with the hard specimen material.

c. Carefully remove the Extended Drill by simultaneously rotating it clockwise while withdrawing it.

Note: To obtain adequate specimen, it may be necessary to repeat the procedure.

OTHER INFORMATION:

Serious incident that has occurred in relation to the device shall be reported in accordance with national laws and regulations. Serious incidents shall be reported to AprioMed AB (legal manufacturer) and to the appropriate regulatory agency such as the competent authority of the Member State (EU/EEA) in which the user/patient is established.

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Rx Only

(EN) Caution: Federal (U.S.A) law restricts this device to sale by or on the order of a physician.

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(EN) Single use. Do not re-use. After use, the product performance and characteristics may degrade, and can cause infection in patient.

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(EN) Do not re-sterilize. Re-sterilization may have impact on product characteristics and performance.

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(EN) Do not use if package is damaged. Device is delivered sterile. Product in damaged package may be contaminated and can cause infection in patient.

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(EN) MR unsafe. Do not use AprioCore plus in Magnetic Resonance Imaging (MRI) equipment. The device is magnetic and exposure to MRI may cause unintended device movement.

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Bonopty® is a registered trademark of AprioMed AB.
Mosaic Surgical Ltd