



IZI Medical Products
From concept to care.



Did you remember to order your MOLDCARE pillows?

MOLDCARE pillows are the perfect compliment to thermoplastic masks for patient immobilization. They provide customized support by molding around a patient's head, neck, and shoulders. When activated with water, the pillow becomes rigid, conforming to the contours of the patient's head and shoulders providing repeatable support for each treatment.

MOLDCARE comes in various sizes to meet your specific needs.

MC7001: Head (standard) 20 cm x 25 cm

MC7002: Head and Shoulder (trapezoid) 58 cm x 65 cm

MC7003: Head (small) 20 cm x 15 cm

MC7004: Head (large) 20 cm x 45 cm

MC7005: Head, Shoulder, Arms 40 cm x 60 cm

MC7006: Head, Shoulder, Arms 60 cm x 60 cm



For Samples or to Place an Order for

MOLDCARE

800.231.1499 or **www.izimed.com**

L089, Rev A