




IZI Medical Products
From concept to care.

Introducing IZI's CushionCast™

IZI Medical Products is proud to introduce the CushionCast™ moldable cushion. CushionCast™ provides a secure custom head support for each patient. Popping the internal pouches releases a resin that hardens the contents of the cushion. The cushion can then be gently formed under your patient's head and neck. It will be a warm pillow that eventually cools to provide secure yet comfortable positioning for highly accurate treatment and reproducibility. The end product is a firm cushion that immobilizes the patient better than the competition and will be dimensionally stable for the entire radiation therapy. CushionCasts™ are designed for use with Silverman head supports or as a support pillow alone, with or without a thermoplastic mask.


CushionCast™ the latest in head immobilization

- Odorless
- Pop internal pouches to release hardening resin
- Once mixed, cushion forms easily and quickly
- Formed cushion is sturdy to fully immobilize
- Larger size increases immobilization



Part Number	Description	Size
CC2335	Moldable Head Cushion	22.5 x 35 cm
CC1520*	Moldable Head Cushion	15.0 x 20 cm

* in development, check with your sales rep for availability



For further information please contact an IZI representative or visit www.izimed.com



IZI Medical Products

From concept to care.

IZI Medical Products Inc. Supplies the following:

ColdCare Pack™
FlexiMarc™ Gold Fiducial Markers
I.Z. PORT® Skin Markers
Indicator® Radiopaque Markers
Indicator® Radiolucent Markers
IZI Gold Fiducial Markers
Klarity® AccuCushion™
Klarity® Oven
Klarity® Thermoplastic Masks
Klarity® Vacuum Bags
Medtronic Probe
MoldCare Pillows
Multi-Modality Fiducial Markers
Multi-Modality Markers
Needle Guides
Quick Spherz® Passive Reflective Markers
ShadowForm® Markers
Snap-On Spherz® Passive Reflective Markers
SoftStrips™ Covers
Threaded Spherz® Passive Reflective Markers
X-MARK™ Gold Fiducial Markers

5 Easter Court, Suite J
Owings Mills, Maryland 21117
T 800.231.1499
T +1.410.594.9403
F +1.410.594.0540