

IZI Coral® Thermoplastic Masks

Stronger, Re-moldable and Reinforced

IZI Medical Products provides our Coral Thermoplastic Masks for strong, fast, comfortable immobilization during radiotherapy of the head and head/neck/shoulders. Special materials prevent sticking, speed the molding process, and enhance reusability. They are available in a variety of thicknesses, perforation types, and apertures for virtually any radiotherapy need.

Product Benefits

- Least shrinkage of all leading thermoplastic masks
- Better working time for a more detailed, better-fitting masks
- 100% non-stick surface when heated
- Stronger than standard RT thermoplastic mask
- T-Lok™ pin design for easy on/off use
- Our masks are reinforced where standard masks are overstretched
- No Latex

Product Specifications

- Ultra-secure reinforcing perforation pattern
- Available in S-Frame and U-Frame

Manufactured in the U.S.



IZI Coral® Thermoplastic Masks Ordering

PART NUMBER	DESCRIPTION	SIZE	THICKNESS	COLOR
MK8007	S-Frame Head Only (reinforced)	7.5"W x 12"H	3.2mm	Coral™
MK8013	U-Frame Head/Neck	7.5"W x 11"H	3.2mm	Coral™
MK8026	Head Only S-Frame	7.5"W x 12"H	3.2mm	Coral™
MK8201	S-Frame, Open Eyes/Mouth (reinforced)	7.5"W x 12.5"H	3.2mm	Coral™
MK8200	S-Frame, Open Eyes/Mouth/Shoulder (reinforced)	21.5"W x 19"H	3.2mm	Coral™