

FlexiMarc™

The New Standard for Gold Fiducial Markers

IZI Medical Products is pleased to introduce the FlexiMarc Fiducial Marker. The first of its kind, this implantable marker produces a unique non-biological node-based projection. The FlexiMarc™ combines 99.9% pure, biocompatible, gold connectors and 3mm long gold nodes. Unlike regular linear or coiled markers, FlexiMarc has nodes that inhibit migration and are easily recognized by the human eye and IGRT software.

Product Benefits

- Uses node technology for better visibility
- Specifically shaped to prevent migration
- Available in 2 or 3 node configurations
- Multiple diameters and needle gauges to meet your needs
- Visible in TRUS, BAT, Portal, MRI, CT, and x-Ray
- 99.9% Gold

Product Specifications

- 0.9mm diameter with 18 GA needle
- 1.2mm diameter with 17 GA needle
- Available in 1cm and 2cm lengths
- Nodes are 3mm in length

Manufactured in the U.S.



FlexiMarc™ Ordering

PART NUMBER	DESCRIPTION
FM09210	FlexiMarc 0.9mm Diameter, 2 Node, 10mm spacing
FM09310	FlexiMarc 0.9mm Diameter, 3 Node, 10mm spacing
FM12210	FlexiMarc 1.2mm Diameter, 2 Node, 10mm spacing
FM12310	FlexiMarc 1.2mm Diameter, 3 Node, 10mm spacing