

I.Z. PORT® Skin Markers

The Most Accurate, Efficient Alignment Markers for Radiation Therapy

Successfully used by thousands of medical centers, our I.Z. PORT® Skin Markers satisfy providers and patients alike. Consisting of a set of markers that outline the treatment field, these markers have a clear backing and split liner that make them easy to apply. Once applied, their thin, accurate lines improve treatment precision. Their ability to remain in place for seven to ten days eliminates the need to re-mark patients and saves time. These markers are also hypoallergenic with no messy ink to enhance patient comfort, and on a clear backing for precise application.

Product Benefits

- ◇ Durability saves times and prevents re-marking
- ◇ Thin, accurate lines for precise treatment
- ◇ Bright blue ink is perfect for all skin types and easy to view
- ◇ Water resistant
- ◇ Sticker-like design is ideal for children

Product Specifications

Lines (IZ1003)

- ◇ 150 lines per package
- ◇ 2.0mm wide x 4.5cm long (0.08" x 1.75")

Isocenters (IZ1004)

- ◇ 300 crosses per package
- ◇ 1.4mm (0.06")

Manufactured in the U.S.



IZ. PORT® Skin Markers Ordering

PART NUMBER	DESCRIPTION
IZ1003	150 IZ Port Lines
IZ1004	300 IZ Port Isocenters