

Indicator® Radiolucent Markers

Cost-effective Markers That Promote A Quick, Effective Diagnosis

IZI Medical Products designed our Indicator® Radiolucent Markers for Mammography with a special pink ink that is radiolucent. Their flat, radiolucent design prevents tissue indentation while facilitating a quick, accurate diagnosis by eliminating artifacts and revealing underlying pathology. The double pull tabs and perforated lines make them easy to apply and remove.

Our markers come in choice of sizes and shapes- including dots for nipples, circles for moles, triangles for palpable masses and lines for scars- to provide the flexibility and accuracy you need. As an added benefit, these markers are disposable and cost effective.

Product Benefits

- Unique pink, radiolucent design
- No artifacts or tissue indentation
- Clear adhesive backing permits exact placement
- Variety of shapes and sizes for all imaging needs
- Flat design
- Double pull tabs
- Disposable

Product Specifications

- 115 dots, circles or triangles per package
- Lines are perforated every 2" for easy use
- Choice of 3 line widths from 0.8mm to 2.5mm

Manufactured in the U.S.



Indicator® Radiolucent Markers Ordering

PART NUMBER	DESCRIPTION
RM2104	115 Dot Markers, 2.5mm
RM2110	66" Line Markers, 1.5mm
RM2111	69" Line Markers, 0.8mm
RM2114	115 Triangle Markers
RM2117	115 Circle Markers