

Indicator® Radiopaque Markers - CT and X-Ray

Precision and Convenience in a Wide Variety of Applications

Our Indicator® Radiopaque Markers can be used in a variety of imaging procedures to identify any point of interest, including:

- Masses
- Scar tissue
- Moles
- Isocenter
- Points of Pain

Users appreciate their ability to provide a clear, accurate reference point that is free of spray artifacts. The flat design prevents tissue indentation and the clear adhesive backing aids in precise marker placement. In addition to enhancing accuracy, they promote procedural efficiency with double pulls tabs that permit easy application and removal. IZI now offers a raised BB marker.

Product Benefits

- Creates precise reference points without artifacts
- Works for CT, X-Ray, fluoroscopy, angiography
- Disposable
- Available in a variety of shapes and sizes for most applicants
- Flat design

Product Specifications

Crosses (TM4001-xx)

- 130 per box, 10 mm and 20 mm sizes, 9.65 mm length

Dots (TM4003-xx)

- 115 per box, 5 sizes from 1.5mm to 4.0mm

Lines (TM4007-xx)

- 62 and 66 each at 2" lengths
- 1.0mm and 1.5mm widths

Arrows (TM4019)

- 115 Markers, 8mm wide

BBs (TM4101-xx)

- 100 Markers with Wings, 1.0mm, 1.5mm, 2.0mm and 2.5mm sizes

BBs (TM4102-xx)

- 100 Markers, 1.0mm, 1.5mm, 2.0mm and 2.5mm sizes

Manufactured in the U.S.


Indicator® Radiopaque Markers- CT and X-Ray Ordering

PART NUMBER	DESCRIPTION	PART NUMBER	DESCRIPTION
TM4001-10	130 Radiopaque Crosses, 10mm	TM4019	115 Radiopaque Arrows
TM4001-20	130 Radiopaque Crosses, 20mm	TM4101-10	100 Marker BB Wings, 1.0mm
TM4003-15	115 Dot Markers, 1.5mm	TM4101-15	100 Marker BB Wings, 1.5mm
TM4003-20	115 Dot Markers, 2.0mm	TM4101-20	100 Marker BB Wings, 2.0mm
TM4003-25	115 Dot Markers, 2.5mm	TM410125	100 Marker BB Wings, 2.5mm
TM4003-30	115 Dot Markers, 3.0mm	TM4102-10	100 Marker BBs, 1.0mm
TM4003-40	115 Dot Markers, 4.0mm	TM4102-15	100 Marker BBs, 1.5mm
TM4007-10	132", 66 Radiopaque Lines, 1.0mm	TM4102-20	100 Marker BBs, 2.0mm
TM4007-15	125", 62 Radiopaque Lines, 1.5mm	TM4102-25	100 Marker BBs, 2.5mm