

Breast Lesion Localization / Jabczenski

Jabczenski Ductogram Cannula

Used for cannulation of the lactiferous duct of the breast

○ The **right-angle design** of the cannula tip enables an **unobstructed view** of the duct for placement and permits the easy **movement and fixation** of the connecting tube away from the area that is being radiographed

Supplied sterile in peel-open packages. Intended for one-time use only.



ORDER INFORMATION

Jabczenski Ductogram Cannula

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	GLOBAL PROD. NO.
DGC-30-5	30G • 5 cm	G04303



TAPERED CANNULA

Not all parts numbers shown may be approved for sale in all regulatory jurisdictions. Parts and configurations are subject to change without notification. Please consult with your Account Representative for details.