

IZI Fiducial Markers

IZI Medical offers a full line of Gold Fiducial Markers and Loose Fiducial Markers. These implantable markers are used with Image Guided Radiation Therapy.

IZI Medical Gold Fiducial Markers and Loose Gold Fiducial Markers are often used in soft tissue such as prostate, lung, liver, and pancreas.

Product Features

- Wide variety of diameters and lengths
- Unique micro-abrasion texturing
- 99.99% pure gold

Product Specifications

- Marker Diameters: 1.2mm, 1.0mm, and 0.8mm
- Marker Lengths: 3mm and 5mm
- Needle Gauges: 17G and 18G
- Needle Lengths: 20cm and 30cm
- Packaging: single pack, three pack, and four pack
- Available loose or in needle
- Tandem markers available

Manufactured in the U.S.



IZI Fiducial Markers Ordering

PART NUMBER	DESCRIPTION
GF0302	2PK, 0.8 x 3.0mm Gold Fiducial, Loose
GF0303	3PK, 0.8 x 3.0mm Gold Fiducial, Loose
GF1001	1PK, 1.2 x 3.0mm Gold Fiducial, 20cm, 17G
GF1003	3PK, 1.2 x 3.0mm Gold Fiducial, 20cm, 17G
GF1004	4PK, 1.2 x 3.0mm Gold Fiducial, 20cm, 17G
GF1321	1PK, 1.0 x 3.0mm Gold Fiducial, 20cm, 18g
GF1323	3PK, 1.0 x 3.0mm Gold Fiducial, 20cm, 18g
GF1324	4PK, 1.0 x 3.0mm Gold Fiducial, 20cm, 18g

PART NUMBER	DESCRIPTION
GF1333	3PK, 1.2 x 3.0mm Gold Fiducial, 30cm, 17G
GF1334	4PK, 1.2 x 3.0mm Gold Fiducial, 30cm, 17G
GF1501	1PK, 1.2 x 5.0mm Gold Fiducial, Loose
GF1503	3PK, 1.2 x 5.0mm Gold Fiducial, Loose
GF1521	1PK, 1.2 x 5.0mm Gold Fiducial, 20cm, 17G
GF1523	3PK, 1.2 x 5.0mm Gold Fiducial, 20cm, 17G
GF1533	3PK, 1.2 x 5.0mm Gold Fiducial, 30cm, 17G
GF1823	3PK, 1.0 x 5.0mm Gold Fiducial, 20cm, 18g