



IZI Medical

Blazer[®]

Vertebral Augmentation System

Curved Cavity Creation
for VCF Treatment

2023 Product Catalog

IZIMed.com



The Blazer[®]-C Vertebral Augmentation System

Minimally Invasive Treatment for Vertebral Compression Fractures

What is The Blazer-C Vertebral Augmentation System?

The Blazer-C Vertebral Augmentation System is a device used in minimally invasive vertebral augmentation procedures that can be performed in a spine specialist's office or in a hospital setting.

The Blazer-C System is indicated for the treatment of pathological compression fractures of the vertebral body that result from osteoporosis, benign lesions, or malignant lesions, by creating channels in the existing spinal bone structure for the flow of bone cement.

How Blazer-C Works

The vertebral augmentation procedure can be done under local or general anesthesia – your doctor will determine which option is best for you. The procedure usually takes less than an hour per fracture treated.


Your doctor will access your VCF via one small (1 cm) incision on your back at the level of the fracture.

- Mechanical channel creation via Nitinol wire
- Predictable cement deposition through cannulated wire
- MIS, Unilateral approach



Blazer-C[®] Case Study

80 YO Female / L2 Fracture / Blazer followed by PMMA Injection



The Blazer-C Vertebral Augmentation System

Mechanical Channel Creation:

The Blazer-C wire allows the physician to control placement of the device and create channels within the vertebral body where desired. Blazer-C can be used to treat pathological vertebral compression fractures including use in sclerotic bone, blastic lesions and chronic fractures.

Predictable Cement Deposition:

Cement flows through the cannulated wire and follows the channels created by Blazer-C providing directional control during cement delivery.

MIS, Unilateral Approach:

The Blazer-C wire is easy to advance and retrieve consistently. The system is designed to provide complete vertebral body coverage via a unilateral approach. It is available in 10G (3mm) outer diameter.



ORDER INFORMATION

Blazer®-C Vertebral Augmentation System

GLOBAL PRODUCT NUMBER

BZK3140	Blazer-C Vertebral Augmentation System with Diamond Trocar Tip Introducer, Radiopaque Bone Cement, Mixing & Injection Kit
ACC3130	Diamond Tip Introducer, 11G • 16 cm
ACC3140	Bevel Tip Introducer, 11G • 16 cm

Not all parts numbers shown may be approved for sale in all regulatory jurisdictions. Parts and configurations are subject to change without notification. Please consult with your Account Representative for details.



Quality Care is our Primary Concern

IZI Medical Products is constantly striving for innovation by continually focusing on developing better products for you, and in the end for the benefit of your patients.

We understand the complex business of healthcare, and we know how important it is to balance quality, with efficiency and affordability.

As a leading innovator, manufacturer, and distributor of quality medical devices, we ensure that all of our products improve patients' outcomes. The result of which, is a better overall experience for you, your team and, most importantly, your patients. We take pride in watching our concepts take shape, transforming and enhancing the way clinical procedures are done. Our innovation and dedication to accuracy is why so many healthcare professionals order our products by name.

We look forward to working together to improve patient care.



IZI Medical

About IZI Medical Products

Based in Owings Mills, Maryland, just outside Baltimore, IZI is a leading innovator, manufacturer, and distributor of quality medical devices that provides the continuity our healthcare partners need to succeed. Over the last three decades, IZI has built a diverse portfolio of products backed by strong intellectual property and currently sells to more than 1,000 domestic customers as well as internationally across 25 countries. IZI's portfolio of high-quality medical consumable devices are used in interventional radiology and oncology, radiation therapy, neuro-spine, and image guided surgery procedures

Our mission is to improve each patient's quality of life by providing our healthcare professional partners a complete portfolio of clinically proven, high value, and innovative medical device solutions inspired by physicians, for physicians. With this mission in mind, our vision is to become a leader in diagnostic and therapeutic minimally invasive solutions to advance healthcare and improve patients' lives.

Contact Us:

Please call 800.231.1499 or +1.410.594.9403 for product samples, or visit IZIMed.com for more information.

We pride ourselves on our responsive service and prompt delivery. When ordered by 4pm Eastern Standard Time (EST), products are shipped same day.

IZI Medical Portfolio also includes products for:

- Biopsy
- Diagnostic Radiology
- Image Guided Surgery
- Radiation Therapy
- Spine Solutions



IZI Medical

5 Easter Court / Suite J
Owings Mills, MD 21117 / USA

P / + 1 800 231 1499
P / + 1 410 594 9403

IZIMed.com