

Better Control. Better Care.

Vertefix® HV. Infused with Insite™ tracking beads.

For vertebroplasty and kyphoplasty procedures, knowing how much and where bone cement is flowing gives you a greater sense of control.

Vertefix HV cement is infused with Insite tracking beads for real-time flow visualization, with high viscosity and long working time, giving you the power to perform and the confidence to deliver better care.

See for yourself the visible difference of Vertefix HV.

What Physicians Are Saying.

“Vertefix HV will help you control your cement and allow you to have excellent cement viscosity immediately after mixing. It has an optimal volume of cement with good working time and viscosity along with great visualization so you can see when the cement is flowing and where it is going. The Insite beads are a great addition and now novel to all of the cements on the market. Bottom line is it’s the best cement I’ve used since I can remember.”

Dr. Douglas Beall
Chief of Services and Interventional Spine Services
Comprehensive Specialty Care, Edmond, Oklahoma


“Improved visualization of cement has always been one of the goals of kyphoplasty and vertebroplasty. If you can visualize the cement better, it makes the procedure safer.”

Dr. Jonathan Morris
Neuroradiologist

Ready To See The Difference?

To learn more about Vertefix HV, and the control that only comes with real-time flow visualization visit:

IZIMed.com/VertefixHV or call us at 800 231 1499



IZI Medical

5 Easter Court / Suite J
Owings Mills, MD 21117 / USA

US / +1 800 231 1499
Int'l / +1 410 594 9403
Info@IZIMed.com

IZIMed.com

Osteo-Site® Products / Cements

Vertefix® HV

Bone Cement Infused with Insite™ Tracking Beads

Used for fixation of pathological fractures of the vertebral body using vertebroplasty and kyphoplasty procedures.

- Two components (powder & liquid) PMMA bone cement
- Cement components are easy to mix, no measurement is required
- Injection time: 18 minutes. Yield: 18 cc of PMMA

Supplied sterile in box. Intended for one-time use.

INDICATIONS

- Painful vertebral compression fractures, benign lesions (hemangioma)
- Malignant lesions (metastatic cancer, myeloma, tumors)
- Treatment of acute & chronic instabilities
- Deformities of the thoracic, lumbar and sacral spine
- Osteoporosis, osteopenia, metastatic disease



Infused with Insite Tracking Beads® for real-time **flow visualization**, with **high viscosity** and long working time for better control.

ORDER INFORMATION

Vertefix® HV

PRODUCT NAME	MIXED VOLUME	CEMENT RADIOPACITY	VISCOSITY	GLOBAL PRODUCT NUMBER
Vertefix HV	18cc	30%	High	G90005

Kit Configurations

PRODUCT NAME	SIZE	DESCRIPTION	GLOBAL PRODUCT NUMBER
VCF Kit	--	DuroJect Injector, MiniMix, and Vertefix HV	G80004
Bone Filler Kit	10G	Osteoforce Injector, Bone Fillers, Syringes, MiniMix, & Vertefix HV	G80007
Bone Filler Kit	11G	Osteoforce Injector, Bone Fillers, Syringes, MiniMix, & Vertefix HV	G80008

Not all parts numbers shown may be approved for sale in all regulatory jurisdictions. Parts and configurations are subject to change without notification. Please consult with your Account Representative for details.