

Osteo-Site® Products / Needles

Murphy Quick Set

Vertebroplasty Set

Used for vertebral body biopsy and infusion

- Includes a connection tube and 10 x 1cc syringes for vertebroplasty injection
- A blunt tip stainless steel obturator is included
- A 30 cm long radiolucent needle holder is included

Supplied sterile in peel-open packages. Intended for one-time use only.

Access and cement injection – all with a **single vertebroplasty kit**



ORDER INFORMATION

Murphy M2™ - MQT

SET ORDER
NUMBER


OUTER CANNULA
GAUGE • LENGTH

GLOBAL
PROD. NO.

DBBN-13-10.0-M2-MQT

13G • 10 cm

G13097



Not all parts numbers shown may be approved for sale in all regulatory jurisdictions. Parts and configurations are subject to change without notification. Please consult with your Account Representative for details.