

Osteo-Site® Products / Needles

Osteo-Site Needle Holder

30 cm needle holder



- Length of holder allows user to fluoroscopically guide needle with reducing radiation exposure
- Radiolucent plastic provides an unobstructed view of needle target and surrounding anatomy
- Compatible with all Osteo-Site Needles

Supplied sterile in peel-open packages. Intended for one-time use only.

Reduces exposure to **radiation** with a secure, non-metallic needle holder.

ORDER INFORMATION

Needle Holder

SET ORDER NUMBER	LENGTH	GLOBAL PROD. NO.
ONH-100	30cm	G13549



Not all parts numbers shown may be approved for sale in all regulatory jurisdictions. Parts and configurations are subject to change without notification. Please consult with your Account Representative for details.