



IZI Medical

izimed.com

Twin Line^{MR}

PLATINUM VISICOIL


Designed for CyberKnife® Image-guided
robotic stereotactic radiosurgery

Now in Platinum for optimal MRI visibility

VISICOIL™

LINEAR FIDUCIAL MARKER

PLATINUM: THE NEW GOLD STANDARD



Benefits of Platinum

ENHANCED MR VISIBILITY

- Better MRI Magnetic Susceptibility
- Has up to 2.5x stronger signal void over gold

“An alternative is to consider magnetic susceptibility when selecting fiducials. This may mean replacing gold with platinum or using a hybrid gold/polymer fiducial.”

MR-Integrated workflows in Radiation Therapy
© Siemens Healthcare GmbH, 2018**

SEE IT!

2 Needles vs. 4 Needles

- Hollow core = minimal CT artifact
- Corrugated surface = excellent ultrasound imaging

TRUST IT!

The Original Linear Marker

- Helical coil design fixates marker in soft tissue
- Stable marker localization from implant through treatment

TREAT IT!

Stranded / Non-Stranded

- Efficient, consistent, and reliable daily treatment set up
- Configured for Cyberknife**

Optimized for CyberKnife®* 2 Markers, 1 Needle



“The Visicoil MR Platinum tracks well with Cyberknife® and the clearer MRI signal void allows for greater confidence with MR/CT fusion for registration. Additionally, the ability to get two of these markers inside one small (19g - 21g) gauge needle is a desirable benefit for our physicians and patients.”

Jun Yang, PhD., DABR Chief Medical Physicist Philadelphia Cyberknife®

CONTACT US TODAY TO LEARN MORE

RadioMed Corporation
Bartlett, Tennessee 38133, USA
North America | T +1 901 432 7202
visicoil@izimed.com | izimed.com

PB-0025 Visicoil-MR-TwinLine_Rev02 | © IZI Medical Products 2020 | All rights reserved | Manufacturer: RadioMed Corporation



*USE OF THE CYBERKNIFE® REGISTERED TRADEMARK BY IZI MEDICAL PRODUCTS IS NOT AN ENDORSEMENT BY INCORPORATED OF THE VISICOIL TWINLINE OR ITS PROMOTION AND ADVERTISING. CYBERKNIFE® IS A TRADEMARK OF ACCURAY INC.
** PUBLISHED BY SIEMENS HEALTHCARE GMBH. ORDER NO. A91MR-1100-95C-7600. PRINTED IN GERMANY. 6085 04180.2