

ShadowForm® Markers

A Better Way to Outline Pelvic Soft Tissue

IZI manufactures disposable ShadowForm® Markers to outline pelvic soft tissue during simulation and CT planning. Patients appreciate the fact that they're soft, smooth, and latex-free while providers like their accuracy and disposability. Safer for high risk patients because there is no cross-contamination. Available for both rectal and vaginal use, these markers are highly flexible and include a T-bar handle that can easily be removed to allow for barium insertion.


Product Benefits

- Available in rectal and vaginal markers
- Materials promote patient comfort and provide ease of use
- Easy barium insertion thanks to removable T-bar handle
- Rectal markers offer excellent measurability
- Prevent Cross Contamination
- Visible on an X-Radiation Modality
- Flexible
- Disposable

Product Specifications

- 10 markers per box
- 1cm intervals on rectal markers
- Vaginal markers (SF5001) 18cm long
- Rectal Markers: (SF5002) 38 cm long, (SF5003) 10 cm long

Manufactured in the U.S.



ShadowForm® Markers Ordering

PART NUMBER	DESCRIPTION	QTY
SF5001	ShadowForm 18cm Vaginal Marker	10 Markers
SF5002	ShadowForm 38cm Rectal Marker	10 Markers
SF5003	ShadowForm 10cm Rectal Marker	10 Markers