




NEW: VISICOIL MR 0.75MM X 0.5CM ON LEFT & 0.75MM X 0.5CM GOLD VISICOIL ON RIGHT



SIEMENS VERIO 3T 70CM OPEN BORE MRI, T2

VISICOILTM MR

|| new ||

Fiducial marker for soft tissue localization

VISICOIL MR in platinum is visible in MRI



IZI Medical[™]

VISICOIL™ MR

Multi-modality fiducial marker

See it!


- Visible under multiple imaging modalities
- Platinum offers enhanced MR visibility compared to gold
- Multiple coil diameters and lengths to suit your needs

Trust it!

- VISICOIL helical design has been available since 2004
- 20 years of clinical data available
- New 21G & 22G needle now available

Ask us for:

- Studies on proton perturbation
- Studies on MRI safety
- Studies on alternative delivery methods





“Gold and Platinum have nearly identical high radiodensities, and are easily seen and tracked by X-ray based imaging. Gold is difficult to see on MRI, whereas Platinum has greater magnetic susceptibility, producing just the right signal void to be readily visible on MRI. Platinum is the ideal fiducial material.”

DR. LEONARD AVRUCH
Department of Medical Imaging / University of Ottawa (Retired)

Available Sizes



Available preloaded in needles or loose



CONTACT US TODAY TO LEARN MORE



IZI Medical™

IZI Medical Products
5 Easter Court / Suite J
Owings Mills, MD 21117 / USA

P / +1 800 231 1499
P / +1 410 594 9403
IZIMed.com