

IZI Gold Fiducial Markers

Precise, Real-Time Localization of Soft Tissue

IZI Medical Products provides all the marking tools you need to better verify localization and reduce your number of fractions. Take our top-selling IZI Gold Fiducial Markers, which enhance accuracy and efficiency during radiation therapy. These implantable markers are used in IMRT, CRT, and IGRT treatment to provide real-time, accurate localization of the treatment volumes. Used in soft tissue cases such as the prostate, liver, and lung, these markers are essential whenever the target may move during and between treatments relative to external skin markings.

Product Benefits

- Ideal size and shape for soft tissue localization
- Unique micro-abrasion texturing reduces markers migration
- Visible in TRUS, BAT, Portal, MRI, CT, and x-Ray
- Superb for image fusion
- 99.9% gold

Product Specifications

- Cylindrical shape
- Sizes: 1.0mm x 3.0mm, 1.2mm x 3.0mm or 1.2mm x 5.0mm
- 17 or 18 gauge needles in 20cm and 30cm lengths
- Available loose or pre-loaded
- **Pre-loaded markers:** Available in packages of 1, 3, and 4
- **Loose markers:** in 3-pack, 0.8mm x 3.0mm
- Tandem dual markers available, 1.2mm x 3.0mm

Manufactured in the U.S.



IZI Gold Fiducial Markers Ordering

PART NUMBER	DESCRIPTION
GF0302	2PK, 0.8 x 3.0mm Gold Fiducial, Loose
GF0303	3PK, 0.8 x 3.0mm Gold Fiducial, Loose
GF1001	1PK, 1.2 x 3.0mm Gold Fiducial, 20cm, 17G
GF1003	3PK, 1.2 x 3.0mm Gold Fiducial, 20cm, 17G
GF1004	4PK, 1.2 x 3.0mm Gold Fiducial, 20cm, 17G
GF1321	1PK, 1.0 x 3.0mm Gold Fiducial, 20cm, 18g
GF1323	3PK, 1.0 x 3.0mm Gold Fiducial, 20cm, 18g
GF1324	4PK, 1.0 x 3.0mm Gold Fiducial, 20cm, 18g

PART NUMBER	DESCRIPTION
GF1333	3PK, 1.2 x 3.0mm Gold Fiducial, 30cm, 17G
GF1334	4PK, 1.2 x 3.0mm Gold Fiducial, 30cm, 17G
GF1501	1PK, 1.2 x 5.0mm Gold Fiducial, Loose
GF1503	3PK, 1.2 x 5.0mm Gold Fiducial, Loose
GF1521	1PK, 1.2 x 5.0mm Gold Fiducial, 20cm, 17G
GF1523	3PK, 1.2 x 5.0mm Gold Fiducial, 20cm, 17G
GF1533	3PK, 1.2 x 5.0mm Gold Fiducial, 30cm, 17G
GF1823	3PK, 1.0 x 5.0mm Gold Fiducial, 20cm, 18g