



IZI Medical

# Quick-Core®

Soft Tissue Biopsy



# Quick-Core<sup>®</sup> Biopsy Needle and Needle Set

Minimally Invasive Soft Tissue Biopsy

## What is Quick-Core<sup>®</sup> Biopsy Needle?

The Quick-Core<sup>®</sup> Biopsy Needle is a device used in minimally invasive soft-tissue biopsy procedures that can be performed in a specialists office or in a hospital setting.

The Quick-Core<sup>®</sup> Biopsy Needle is indicated for soft-tissue biopsy.

## How Quick-Core<sup>®</sup> Biopsy Needle and Needle Set Works

The soft tissue biopsy procedure can be done under local anesthetic. The procedure usually takes less than an hour to complete.

The Quick-Core<sup>®</sup> set includes a stainless steel trocar.

## Why use Quick-Core<sup>®</sup> Biopsy Needle versus Competitors?

- Smart design reduces moving components allowing multiple uses without breakage
- Sharp bevel tip with two additional facets facilitates easy placement into tissue
- Rigid spring reduces vibration which allows better cutting of specimen preserving tissue architecture
- Small, ergonomic and lightweight handle limits needle movement
- No device recalls and a legacy of quality



Quick-Core® / Needle


# Quick-Core® Biopsy Needle

## Soft-tissue Biopsy Needle

Used for soft-tissue biopsy

- The lightweight and compact Quick-Core® Biopsy Needle **facilitates one-handed operation** during procedures
- The **precise, rapid-firing mechanism** provides an intact core tissue sample
- The beveled stylet **permits easy penetration** into the lesion with less trauma to surrounding tissues and organs
- The etched tip **enhances visualization** under ultrasound


Supplied sterile in peel-open packages. Intended for one-time use only.



ORDER INFORMATION

**Quick-Core® Biopsy Needle**

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	GLOBAL PROD. NO.
<b>10 MM THROW LENGTH</b>		
QC-14-6.0-10T	14G • 6 cm	G28104
QC-14-9.0-10T	14G • 9 cm	G28105
QC-14-11.0-10T	14G • 11 cm	G28127
QC-14-15.0-10T	14G • 15 cm	G28115
QC-14-20.0-10T	14G • 20 cm	G28116
QC-16-6.0-10T	16G • 6 cm	G28280
QC-16-9.0-10T	16G • 9 cm	G28282
QC-16-11.0-10T	16G • 11 cm	G28129
QC-16-15.0-10T	16G • 15 cm	G28110
QC-16-20.0-10T	16G • 20 cm	G28118
QC-18-6.0-10T	18G • 6 cm	G28283
QC-18-9.0-10T	18G • 9 cm	G27009
QC-18-11.0-10T	18G • 11 cm	G28131
QC-18-15.0-10T	18G • 15 cm	G07819
QC-18-20.0-10T	18G • 20 cm	G28120
QC-20-6.0-10T	20G • 6 cm	G27012
QC-20-9.0-10T	20G • 9 cm	G28126
QC-20-11.0-10T	20G • 11 cm	G28133
QC-20-15.0-10T	20G • 15 cm	G07822
QC-20-20.0-10T	20G • 20 cm	G42020
<b>20 MM THROW LENGTH</b>		
QC-14-6.0-20T	14G • 6 cm	G28100
QC-14-9.0-20T	14G • 9 cm	G27006
QC-14-11.0-20T	14G • 11 cm	G28107
QC-14-15.0-20T	14G • 15 cm	G07824
QC-14-20.0-20T	14G • 20 cm	G28102
QC-16-6.0-20T	16G • 6 cm	G28281
QC-16-9.0-20T	16G • 9 cm	G27007
QC-16-11.0-20T	16G • 11 cm	G28117
QC-16-15.0-20T	16G • 15 cm	G07794
QC-16-20.0-20T	16G • 20 cm	G07827
QC-18-6.0-20T	18G • 6 cm	G27008
QC-18-9.0-20T	18G • 9 cm	G27010
QC-18-11.0-20T	18G • 11 cm	G28119
QC-18-15.0-20T	18G • 15 cm	G07778
QC-18-20.0-20T	18G • 20 cm	G07821
QC-20-6.0-20T	20G • 6 cm	G28124
QC-20-9.0-20T	20G • 9 cm	G27015
QC-20-11.0-20T	20G • 11 cm	G28122
QC-20-15.0-20T	20G • 15 cm	G07780
QC-20-20.0-20T	20G • 20 cm	G42021



THROW LENGTH

Not all parts numbers shown may be approved for sale in all regulatory jurisdictions. Parts and configurations are subject to change without notification. Please consult with your Account Representative for details.

Quick-Core® / Set


## Quick-Core® Biopsy Needle Set

### Soft-tissue Biopsy Needle Set

Used for soft-tissue biopsy

- Coaxial outer needle permits multiple biopsies with the Quick-Core® needle through a **single access site**
- The lightweight and compact Quick-Core® Biopsy Needle facilitates **one-handed operation** during procedures
- The **precise, rapid-firing mechanism** provides an intact core tissue sample
- The EchoTip® surface of the outer coaxial needle and the etched tip of the cutting cannula **enhance visibility** under ultrasound


Supplied sterile in peel-open packages. Intended for one-time use only.



ORDER INFORMATION

**Quick-Core® Biopsy Needle Set**

SET ORDER NUMBER	BIOPSY NEEDLE GAUGE • LENGTH	COAXIAL NEEDLE GAUGE • LENGTH	GLOBAL PROD. NO.
<b>10 MM THROW LENGTH</b>			
QCS-14-9.0-10T	14G • 9 cm	11G • 6.3 cm	G28106
QCS-14-11.0-10T	14G • 11 cm	11G • 8.4 cm	G28128
QCS-14-15.0-10T	14G • 15 cm	11G • 12.3 cm	G28108
QCS-14-20.0-10T	14G • 20 cm	11G • 17.3 cm	G28109
QCS-16-9.0-10T	16G • 9 cm	14G XTW • 6.5 cm	G28111
QCS-16-11.0-10T	16G • 11 cm	14G XTW • 8.6 cm	G28130
QCS-16-15.0-10T	16G • 15 cm	14G XTW • 12.5 cm	G28112
QCS-16-20.0-10T	16G • 20 cm	14G XTW • 17.5 cm	G28114
QCS-18-9.0-10T	18G • 9 cm	16G TW • 6 cm	G27016
QCS-18-11.0-10T	18G • 11 cm	16G TW • 8 cm	G28132
QCS-18-15.0-10T	18G • 15 cm	16G TW • 12 cm	G08728
QCS-18-20.0-10T	18G • 20 cm	16G TW • 16.9 cm	G28121
QCS-20-9.0-10T	20G • 9 cm	18.5G XTW • 6.5 cm	G27018
QCS-20-11.0-10T	20G • 11 cm	18.5G XTW • 8.5 cm	G28134
QCS-20-15.0-10T	20G • 15 cm	18.5G XTW • 12.5 cm	G08727
QCS-20-20.0-10T	20G • 20 cm	18.5G XTW • 17.4 cm	G28123
<b>20 MM THROW LENGTH</b>			
QCS-14-9.0-20T	14G • 9 cm	11G • 5.3 cm	G28290
QCS-14-11.0-20T	14G • 11 cm	11G • 8.4 cm	G51411
QCS-14-15.0-20T	14G • 15 cm	11G • 11.3 cm	G28101
QCS-14-20.0-20T	14G • 20 cm	11G • 16.3 cm	G28103
QCS-16-9.0-20T	16G • 9 cm	14G XTW • 5.5 cm	G28292
QCS-16-11.0-20T	16G • 11 cm	14G XTW • 8.6 cm	G51611
QCS-16-15.0-20T	16G • 15 cm	14G XTW • 11.5 cm	G10149
QCS-16-20.0-20T	16G • 20 cm	14G XTW • 16.5 cm	G28113
QCS-18-9.0-20T	18G • 9 cm	16G TW • 5 cm	G27017
QCS-18-11.0-20T	18G • 11 cm	16G TW • 8 cm	G51811
QCS-18-15.0-20T	18G • 15 cm	16G TW • 11 cm	G08790
QCS-18-20.0-20T	18G • 20 cm	16G TW • 16 cm	G08789
QCS-20-9.0-20T	20G • 9 cm	18.5G XTW • 5.5 cm	G27019
QCS-20-11.0-20T	20G • 11 cm	18.5G XTW • 8.5 cm	G52011
QCS-20-15.0-20T	20G • 15 cm	18.5G XTW • 11.5 cm	G08788
QCS-20-20.0-20T	20G • 20 cm	18.5G XTW • 16.5 cm	G08787



XTW- EXTRA THIN WALL  
TW- THIN WALL

Not all parts numbers shown may be approved for sale in all regulatory jurisdictions. Parts and configurations are subject to change without notification. Please consult with your Account Representative for details.



## Quality Care is our Primary Concern

IZI Medical Products is constantly striving for innovation by continually focusing on developing better products for you, and in the end for the benefit of your patients.

**We understand the complex business of healthcare, and we know how important it is to balance quality, with efficiency and affordability.**

As a leading innovator, manufacturer, and distributor of quality medical devices, we ensure that all of our products meet or exceed U.S. and international regulatory and quality standards. The result of which, is a better overall experience for you, your team and, most importantly, your patients. We take pride in watching our concepts take shape, transforming and enhancing the way clinical procedures are done. Our innovation and dedication to accuracy is why so many healthcare professionals order our products by name.

We look forward to working together to improve patient care.





IZI Medical

## About IZI Medical Products

Based in Owings Mills, Maryland, just outside Baltimore, IZI is a leading innovator, manufacturer, and distributor of quality medical devices that provides the continuity our healthcare partners need to succeed. Over the last three decades, IZI has built a diverse portfolio of products backed by strong intellectual property and currently sells to more than 1,000 domestic customers as well as internationally across 25 countries. IZI's portfolio of high-quality medical consumable devices are used in interventional radiology and oncology, radiation therapy, neuro-spine, and image guided surgery procedures

Our mission is to improve each patient's quality of life by providing our healthcare professional partners a complete portfolio of clinically proven, high value, and innovative medical device solutions inspired by physicians, for physicians. With this mission in mind, our vision is to become a leader in diagnostic and therapeutic minimally invasive solutions to advance healthcare and improve patients' lives.


Contact Us:

Please call 800.231.1499 or +1.410.594.9403 for product samples, or visit [IZIMed.com](http://IZIMed.com) for more information.

We pride ourselves on our responsive service and prompt delivery. When ordered by 4pm Eastern Standard Time (EST), products are shipped same day.

IZI Medical Portfolio also includes products for:

- Biopsy
- Diagnostic Radiology
- Image Guided Surgery
- Radiation Therapy
- Spine Solutions
- Breast Lesion Localization





**IZI Medical**

**UNITED STATES**

**IZI MEDICAL PRODUCTS, LLC**

Attn: Customer Service Department  
5 Easter Court, Suite J  
Owings Mills, Maryland 21117  
United States of America  
Phone: +1 (410) 594-9403  
Fax: +1 (410) 594-0540  
Email: [orders@izimed.com](mailto:orders@izimed.com)

---

**ITALY**

**MED-ITALIA BIOMEDICA, S.r.L**

Via Piero Gobetti, 56r.  
16145 Genova  
Italy  
P / + 39 010 3080301  
F / + 39 010 3080311  
[customerservice@meditalia.net](mailto:customerservice@meditalia.net)

---

**BELGIUM & LUXEMBURG**

**RMS MEDICAL DEVICES**

Omer De Vidtslaan 50  
1760 Roosdaal  
Belgium  
P / + 32 (0) 54 30 01 10  
F / + 32 (0) 54 34 33 40  
[sales@rmsmedicaldevices.com](mailto:sales@rmsmedicaldevices.com)

---

**SWITZERLAND, AUSTRIA & GREECE**

**STÖCKLI MEDICAL AG**

Länggasse 4  
6208 Oberkirch  
Switzerland  
P / + 41 41 925 66 55  
F / + 41 41 925 66 77  
[Info@stoeklimedical.ch](mailto:Info@stoeklimedical.ch)

---

**SPAIN**

**PRIM S.A.**

Calle Yolanda González, número 15  
Polígono Industrial número 1  
28938 - Móstoles, Madrid  
P / + 34 91 334 24 00  
F / + 34 91 334 24 94  
[atencion.cliente@prim.es](mailto:atencion.cliente@prim.es)

**FRANCE**

**TSC, THE SURGICAL  
COMPANY**

Attn: Customer Service Department  
26 Rue de la Montée  
68720 Flaxlanden  
France  
P / + 33 03 89 06 14 44  
F / + 33 03 89 06 41 99  
[info@tscfrance.com](mailto:info@tscfrance.com)

---

**UNITED KINGDOM**

**UK MEDICAL**

Albreda House  
Lydgate Lane  
Sheffield  
S10 5FH  
United Kingdom  
P / + 44 (0) 114 268 8880  
F / + 44 (0) 114 268 8881  
[hello@ukmedical.com](mailto:hello@ukmedical.com)

---

**NETHERLANDS**

**RMS MEDICAL DEVICES**

Flight Forum 40 Ground Floor  
5657 DB Eindhoven  
The Netherlands  
P / + 31 61 23 51 821  
F / + 32 (0) 54 34 33 40  
[sales@rmsmedicaldevices.com](mailto:sales@rmsmedicaldevices.com)

**JAPAN**

**SUGAN CO., LTD.**

1-12-10 Higashi Tenma  
Kita-ku, Osaka 530-0044  
Japan  
P / + 06 6354 6181  
F / + 06 6354 8007

---

**AUSTRALIA**

**POINT BLANK MEDICAL**

F12, 1-15 Barr Street  
Balmain  
New South Wales 2041  
Australia  
P / + 61 (0) 456 438 409  
F / + 61 (02) 8088 6595  
[enquiries@pointblankmedical.com](mailto:enquiries@pointblankmedical.com)

---

**ALL OTHER EUROPEAN, MIDDLE EASTERN  
& AFRICAN COUNTRIES**

**ILUMARK GmbH**

Hohenlindnerstr. 11c  
85622 Feldkirchen  
Germany  
P / + 49 89 907795980  
F / + 49 89 907795989  
[info@ilumark.com](mailto:info@ilumark.com)

---

**NORDIC COUNTRIES**

**NORDIC DISTRIBUTOR**

Vingmed-group  
[www.vingmed.com](http://www.vingmed.com)



**IZI Medical**

5 Easter Court / Suite J  
Owings Mills, MD 21117 / USA

P / + 1 800 231 1499  
P / + 1 410 594 9403

[IZIMed.com](http://IZIMed.com)