



Breast Lesion Localization / Müller-Schimpfle Coil

Müller-Schimpfle Breast Localization Coil

Compatible with Ultrasound, CT and MRI

Used for preoperative marking of non palpable breast lesions including axillary lesions

- The EchoTip® needle allows for **increased visibility** under ultrasound
- The coil acts as a **marker** to indicate the location of the lesion
- The device is MR conditional up to **3 tesla**

Supplied sterile in peel-open packages. Intended for one-time use only.

ORDER INFORMATION

Breast Localization Coil

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	COIL DIAMETER	GLOBAL PROD. NO.
IDLC-19.5-9.0-4-U-DM	19.5G • 9 cm	4mm	G20782

