

Needle Guides

Enhancing Accuracy and Speed in CT-Guided Procedures

IZI Medical Products' specially designed Needle Guides reduce procedure time and increase precision in CT-guided pain injections and biopsies. These flat, self-adhesive markers aid in determining needle angle and entry site by providing a clear, accurate reference point with no spray artifacts. IZI Needle Guides are used by neurosurgeons, neuroradiologists and interventional radiologists at leading health institutions around the world.

Product Benefits

- O Flat, disposable markers do not indent tissue
- O Clear adhesive backing permits visualization of exact marker placement
- Easy-to-read pattern helps identify accurate needle angle and entry site
- Windows allow needle passage
- Wider first line spacing aids orientation
- Thicker center line facilitates counting
- Flexible to contour to the body
- O Choice of two sizes provides options for various locations and body types

Product Specifications

Large (NG4010-1)

- O Ideal for broader areas such as the back
- O Grid measures 3"H x 4"W
- Pattern: series of 7 lines creates8 access windows
- 15 grids per box

Small (NG4010-2)


- O Ideal for smaller areas such as the neck
- Orid measures 1.5"H x 4"W
- Pattern: series of 7 lines creates8 access windows
- 15 grids per box

For a video, please visit: IZIMed.com/needleguide

Manufactured in the U.S.

13 2019 Product Guide





Needle Guides Ordering

PART NUMBER	DESCRIPTION	QUANTITY
NG4010-1	Large Needle Guide	15 Grids / Box
NG4010-2	Small Needle Guide	15 Grids / Box