

## Disposable Passive Array - Pedicle Access Needle

### Image Guided Surgery

IZI Medical Products' Disposable Passive Array - Pedicle Access Needle has been specifically designed for use with Medtronic's\* image guided surgery (IGS) systems, allowing accurate access through the pedicle. IZI's Pedicle Access Needles are sterile single-use instruments providing increased safety, time savings, and reduced costs.

### Security

All components of the navigated pedicle access kit exceed highest quality standards. The stainless-steel cannula is rigidly connected to a stainless steel Luer-lock connector. The flattened area on the handle allows easier mallet access, so impact force is directed straight into the high precision facet cut tip of the needle. With an easy rotation movement the needle can be removed from the cannula, immediately after unlocking.

### Quality

Our Navigation Markers enable unlimited and precise visibility of the Pedicle Access Needle by the IGS-System.

- Precise alignment of pedicle access
- Intuitive and swift locking of cannula
- Extra-long, rigid cannula
- Ergonomic, lightweight handle
- Flat surface for easy mallet access

### Range of Application

The extended cannula length (160mm) allows the use of the Pedicle Access Needle even with larger patients.

### Product Specifications

- 1 Pedicle Access Needle with handle and Navigation array includes 4 Navigation Markers
- 2 cannulas (ID 2.1 mm, length 208 mm)
- 1 removable handle

### Control

Cannula locking is ergonomically designed using a visual and haptic lock/unlock control.

\*Medtronic\* and StealthStation\* are registered trademarks of Medtronic, Inc. Medtronic does not endorse nor is affiliated with IZI Medical Products.

Manufactured in the U.S.



## Disposable Passive Array - Pedicle Access Needle Ordering

PART NUMBER	DESCRIPTION
-------------	-------------

MP6001	Contents per unit:
--------	--------------------

	1 Pedicle Access Needle with Navigation array including 4 Navigation markers
--	--

	2 cannulas (ID 2.1mm, length 208mm)
--	-------------------------------------

	1 removable handle
--	--------------------