



IZI Medical™

VISICOIL™

FIDUCIAL MARKER

Visicoil.com


IZIMed.com



Products for Peak Radiation Therapy Performance.

Radiation oncologists, therapists, oncology nurses, and technologists like you understand the value of time and accuracy. That's why **IZI Medical Products®** ensures that our tools reduce procedural time and help you target malignant tissue, creating a better overall experience for you, your team, and your patients.

We take pride in watching our concepts become products that transform and enhance your ability to perform the highest quality simulation and radiotherapy treatments. From simulation through treatment, our products allow you to provide optimal patient stabilization and comfort while providing all the marking tools you need to better verify localization and reduce your number of fractions. Our innovation and dedication to accuracy are why so many healthcare professionals order our radiotherapy products by name.




Quality Care is our Primary Concern

IZI Medical Products is constantly striving for innovation by continually focusing on developing better products for you, and in the end for the benefit of your patients.

We understand the complex business of healthcare, and we know how important it is to balance quality, with efficiency and affordability.

As a leading innovator, manufacturer, and distributor of quality medical devices, we ensure that all of our products reduce procedural time and enhance precision. The result of which, is a better overall experience for you, your team and, most importantly, your patients. We take pride in watching our concepts take shape, transforming and enhancing the way clinical procedures are done. Our innovation and dedication to accuracy is why so many healthcare professionals order our products by name.



Visicoil™ Fiducial Markers

Fiducial Markers for Soft Tissue Tumor Localization

The VISICOIL™ / VISICOIL^{MR} are available preloaded in various needle gauges with multiple marker sizes and lengths. This helps you choose the right implant method for your patient. The correct marker will optimally reduce common image distortions.

Product Benefits

- Wide range of needle sizes for physician preference
- Centimeter marks for depth placement
- Echogenic needle tip
- Pre-Loaded, packaged sterile and ready for use
- Platinum has Improved MRI visibility compared to gold

Product Specifications

- Platinum Marker Diameter: 0.35mm, 0.50mm, 0.75mm, 1.00mm
- Platinum Marker Length: 5mm, 10mm
- Gold Marker Diameter: 0.35mm, 0.50mm, 0.75mm, 1.10mm
- Gold Marker Length: 5mm, 10mm, 20mm, 30mm
- Needle Gauge: 17g, 18g, 19g, 20g
- Needle Length: 12cm, 20cm, 30cm

ORDER INFORMATION

Platinum

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
PRE-LOADED		
MR-035-005-PL	19G • 20cm	0.35mm • 5mm
MR-035-010-PL	19G • 20cm	0.35mm • 10mm
MR-035-020-PL	19G • 20cm	0.35mm • 20mm
MR-050-005-1920	19G • 20cm	0.50mm • 5mm
MR-050-010-1920	19G • 20cm	0.50mm • 10mm
MR-075-005-PL	19G • 20cm	0.75mm • 5mm
MR-075-010-PL	19G • 20cm	0.75mm • 10mm
MR-100-005-PL	18G • 20cm	1.00mm • 5mm
MR-100-010-PL	18G • 20cm	1.00mm • 10mm

ORDER INFORMATION

Gold

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
PRE-LOADED		
VC-035-005-PL	19G • 20 cm	0.35mm • 5mm
VC-035-010-PL	19G • 20 cm	0.35mm • 10mm
VC-035-020-PL	19G • 20cm	0.35mm • 20mm
VC-050-005-PL	18G • 20cm	0.50mm • 5mm
VC-050-010-PL	18G • 20cm	0.50mm • 10mm

ORDER INFORMATION

Gold

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
PRE-LOADED		
VC-075-005-PL	18G • 20cm	0.75mm • 5mm
VC-075-010-PL	18G • 20cm	0.75mm • 10mm
VC-075-020-PL	18G • 20cm	0.75mm • 20mm
VC-075-030-PL	18G • 20cm	0.75mm • 30mm
VC-110-005-PL	17G • 20cm	1.10mm • 5mm
VC-110-010-PL	17G • 20cm	1.10mm • 10mm
VC-110-020-PL	17G • 20cm	1.10mm • 20mm
VC-050-005-PLSL	19G • 20cm	0.50mm • 5mm
VC-050-010-PLSL	19G • 20cm	0.50mm • 10mm
VC-050-020-PLSL	19G • 20cm	0.50mm • 20mm
VC-050-005-PL30	18G • 30cm	0.50mm • 5mm
VC-050-010-PL30	18G • 30cm	0.50mm • 10mm
VC-075-005-PL12	18G • 12cm	0.75mm • 5mm
VC-075-010-PL12	18G • 12cm	0.75mm • 10mm
VC-075-005-PL30	18G • 30cm	0.75mm • 5mm
VC-075-010-PL30	18G • 30cm	0.75mm • 10mm
VC-075-020-PL30	18G • 30cm	0.75mm • 20mm
VC-110-010-PL30	17G • 30cm	1.10mm • 10mm
VC-110-020-PL30	17G • 30cm	1.10mm • 20mm



Visicoil™ Fiducial Markers

Fiducial Markers for Soft Tissue Tumor Localization

VISICOIL™ / VISICOIL^{MR} 0.50mm and 0.35mm markers are available pre-loaded in 21 and 22 gauge needles. Commonly used for implants in Interventional Radiology or anywhere thinner gauge needles are preferred.

Product Benefits

- 50% reduction in gauge size vs. traditional gold fiducials
- Shorter 8cm and 15cm needles available
- Single bevel needle
- Centimeter marks for depth placement

Product Specifications

- Platinum Marker Diameter: 0.35mm, 0.50mm
- Platinum Marker Length: 5mm, 10mm, 20mm
- Gold Marker Diameter: 0.35mm, 0.50mm
- Gold Marker Length: 5mm, 10mm, 20mm
- Needle Gauge: 21g, 22g
- Needle Length: 8cm, 15cm, 20cm



ORDER INFORMATION

Platinum

SET ORDER NUMBER	NEEDLE GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
INTERVENTIONAL RADIOLOGY		
SNMR-035-005-2208	22G • 8cm	0.35mm • 5mm
SNMR-035-005-2215	22G • 15cm	0.35mm • 5mm
SNMR-035-005-2220	22G • 20cm	0.35mm • 5mm
SNMR-035-010-2208	22G • 8cm	0.35mm • 10mm
SNMR-035-010-2215	22G • 15cm	0.35mm • 10mm
SNMR-035-010-2220	22G • 20cm	0.35mm • 10mm
SNMR-035-020-2215	22G • 15cm	0.35mm • 20cm
SNMR-035-020-2220	22G • 20cm	0.35mm • 20mm
SNMR-050-005-2108	21G • 8cm	0.50mm • 5mm
SNMR-050-005-2115	21G • 15cm	0.50mm • 5mm
SNMR-050-005-2120	21G • 20cm	0.50mm • 5mm
SNMR-050-010-2108	21G • 8cm	0.50mm • 10mm
SNMR-050-010-2115	21G • 15cm	0.50mm • 10mm
SNMR-050-010-2120	21G • 20cm	0.50mm • 10mm

ORDER INFORMATION

Gold

SET ORDER NUMBER	NEEDLE GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
INTERVENTIONAL RADIOLOGY		
SN-035-005-PL08	22G • 8cm	0.35mm • 5mm
SN-035-005-PL15	22G • 15cm	0.35mm • 5mm
SN-035-010-PL08	22G • 8cm	0.35mm • 10mm
SN-035-010-PL15	22G • 15cm	0.35mm • 10mm
SN-035-020-PL15	22G • 15cm	0.35mm • 20mm
SN-050-005-PL08	21G • 8cm	0.50mm • 5mm
SN-050-005-PL15	21G • 15cm	0.50mm • 5mm
SN-050-010-PL08	21G • 8cm	0.50mm • 10mm
SN-050-010-PL15	21G • 15cm	0.50mm • 10mm
SN-050-020-PL15	21G • 15cm	0.50mm • 20mm



Visicoil™ TwinLine Fiducial Markers

Linear Fiducial Markers for Soft Tissue Tumor Localization

VISICOIL™ TwinLine and VISICOIL^{MR} TwinLine are designed for image-guided robotic stereotactic radiosurgery. TwinLine stranded and non-stranded contain 2 VISICOIL™ markers with a specific length spacer, pre-loaded, sterile and ready to use.

Product Benefits

- Designed for Cyberknife* image-guided robotic stereotactic radiosurgery
- Two markers, one needle
- Consistent marker spacing
- Stranded / Non-stranded
- Platinum offers improved MRI visibility compared to gold

Product Specifications

- Marker Diameter: 0.50mm, 0.75mm
- Marker Length: 5mm
- Spacer Non-Stranded: 15mm, 20mm
- Spacer Stranded: 15mm, 20mm
- Needle Gauge: 18g, 19g, 20g, 21g
- Needle Length: 12cm, 20cm

Use of the Cyberknife registered trademark by IZI Medical Products is not an endorsement by incorporated of the Visicoil TwinLine or its promotion and advertising. Cyberknife* is a trademark of Accury Inc..



ORDER INFORMATION



Platinum

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH	SPACER LENGTH	OUTER SLEEVE
TWINLINE				
MRTL-050-41520	19G • 20cm	0.50mm • 5mm	15mm	
MRTL-050-41521	21G • 20cm	0.50mm • 5mm	15mm	
PTTL-050-11021	21G • 20cm	0.50mm • 5mm	10mm	
MRTLS-075-51518	18G • 20cm	0.75mm • 5mm	15mm	Yes

ORDER INFORMATION

Gold

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH	SPACER LENGTH	OUTER SLEEVE
TWINLINE				
TL-050-21520	20G • 20cm	0.50mm • 5mm	15mm	
TL-075-21518	18G • 20cm	0.75mm • 5mm	15mm	
TL-050-62020	20G • 20cm	0.50mm • 5mm	20mm	
TLS-075-31518	20G • 20cm	0.75mm • 5mm	15mm	Yes
TLS-075-31518-12	18G • 12cm	0.75mm • 5mm	15mm	Yes



Visicoil™ Fiducial Markers

Fiducial Markers for Soft Tissue Tumor Localization

IZI Medical Products offers VISICOIL™ for use with a variety of image guided radiotherapy treatment plans. VISICOIL™/VISICOIL^{MR} are unique helically wound 99.99% pure gold or 99.95% pure platinum fiducial markers. VISICOIL™ has been available since 2004 and has been used in hundreds of thousands of patients worldwide. Platinum has a higher magnetic susceptibility.

Product Benefits

- Corrugated Surface
- Helically wound with a hollow core
- Elongated design may allow for fewer markers being required for IGRT
- Platinum = Improved MRI visibility

Product Specifications

- Platinum Marker Diameter: 0.35mm, 0.50mm, 0.75mm, 1.00mm
- Platinum Marker Length: 5mm, 10mm
- Gold Marker Diameter: 0.35mm, 0.50mm, 0.75mm, 1.10mm
- Gold Marker Length: 5mm, 10mm, 20mm, 30mm





ORDER INFORMATION

Platinum

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
LOOSE		
MRFP-035-005		0.35mm • 5mm
MRFP-035-010		0.35mm • 10mm
MRFP-050-005		0.50mm • 5mm
MRFP-050-010		0.50mm • 10mm
MRFP-075-005		0.75mm • 5mm
MRFP-075-010		0.75mm • 5mm
MRFP-100-005		1.00mm • 5mm
MRFP-100-010		1.00mm • 10mm

ORDER INFORMATION

Gold

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
LOOSE		
FP-03505		0.35mm • 5mm
FP-0009		0.35mm • 10mm
FP-05005		0.50mm • 5mm
FP-05010		0.50mm • 10mm
FP-0755		0.75mm • 5mm
FP-0015		0.75mm • 10mm
FP-11005		1.10mm • 5mm
FP-0057		1.10mm • 10mm

About IZI Medical Products

Based in Owings Mills, Maryland, just outside Baltimore, IZI is a leading innovator, manufacturer, and distributor of quality medical devices that provides the continuity our healthcare partners need to succeed. Over the last three decades, IZI has built a diverse portfolio of products backed by strong intellectual property and currently sells to more than 1,000 domestic customers as well as internationally across 25 countries. IZI's portfolio of high-quality medical consumable devices are used in interventional radiology and oncology, radiation therapy, neuro-spine, and image guided surgery procedures

Our mission is to improve each patient's quality of life by providing our healthcare professional partners a complete portfolio of clinically proven, high value, and innovative medical device solutions inspired by physicians, for physicians. With this mission in mind, our vision is to become a leader in diagnostic and therapeutic minimally invasive solutions to advance healthcare and improve patients' lives.

Please call for product samples, or visit izimed.com for more information.

- Radiation Therapy
- Diagnostic Radiology
- Image Guided Surgery
- Biopsy Localization
- Spine Solutions
- Breast Lesion Localization
- Biopsy
- Diagnostic Radiology



IZI Medical Products
5 Easter Court / Suite J
Owings Mills, MD 21117 / USA

P / + 1 800 231 1499
P / + 1 410 594 9403
IZIMed.com